



**PHILIPS**

Radiation Oncology

MRI

## What's the right MRI field strength for radiation therapy planning?

MR scanners are becoming an integral part of radiation oncology workflows as hospitals experience the benefits of using MRI for radiation therapy planning. But do you need a 1.5T or 3.0T MRI system?

To help you answer this question, we turned to Marielle Philippens, MD, PhD, who has experience using both. As a medical physicist in the Department of Radiotherapy at the University Medical Center in Utrecht (UMCU), the Netherlands, Dr. Philippens has worked with several Philips Ingenia MR-RT (1.5T and 3.0T) systems since 2013.

# Do you need a 1.5T or 3.0T MRI system?

## Site-specific patient population

“The high 3.0T field strength provides higher SNR and spatial resolution than 1.5T. However, at 1.5T distortions tend to be fewer than at 3.0T. The choice between 1.5T and 3.0T can depend on the anatomies that will need to be scanned and the desired balance between the required resolution and the ability to interpret distortions. For instance, the head-neck area is more sensitive to distortions than the prostate.”

## MRI knowledge and experience of the staff

“3.0T MRI requires a higher staff expertise than 1.5T. So, for an oncology department starting out in using MRI for radiation therapy planning, choosing 1.5T can help the team become accustomed to the new workflow.”

## Treatment planning only or broader use

“If treatment simulation is the main use, I recommend to consider starting with 1.5T, as it is easier for reducing distortions, although resolution and contrast for brain and pelvic tumors are not as good as with 3.0T. If MRI is also used for response monitoring and functional imaging of the tumor, I think 3.0T should be the first choice, although it has to be considered that thoracic imaging (esophagus and lung) are not feasible on 3.0T.”

## Combination with 1.5T MR-linac

“For simulation of a treatment at a 1.5T MR-linac, 1.5T is the most obvious choice as contrast on this scanner is similar to the contrast on the MR-linac. In this way, the same imaging parameters can be used on the MR simulator and on the MR-linac. It also helps to estimate susceptibility effects on the geometric accuracy of the images. A 3.0T scanner can be used as well to simulate MR-linac treatments without major disadvantages.”

## When both 1.5T and 3.0T are available

“If a hospital has access to both 1.5T and 3.0T, I'd recommend doing MRI for prostate and rectal cancer planning on 3.0T, and asking the radiologist for a preference on other applications, because the choice depends on the patient and the disease,” says Dr. Philippens. “It also depends on the knowledge of the staff and the medical physicist, who need to adjust the ExamCards and sequences between the two field strengths.”

This information is extracted from the FieldStrength publication [“Approaches for including MRI in radiation therapy planning”](#).



**Marielle Philippens,  
MD, PhD**

Marielle Philippens, MD, PhD, has been a radiobiologist and medical physicist with focus on MRI physics in the department of Radiation Oncology at University Medical Center Utrecht since 2008. Her research interest is functional imaging for oncology, with particular focus on diffusion weighted MR imaging. Her research focus areas are head and neck cancer, rectal cancer and breast cancer

### Differences between field strengths

This table summarizes the view of the radiation therapy department at UMC Utrecht on differences between 1.5T or 3.0T field strengths for use in radiation therapy applications.

Acquisition techniques	1.5T	3.0T
T2W TSE	Voxel size: 0.7 x 0.7 x 3 mm <sup>3</sup>	Voxel size: 0.45 x 0.45 x 3 mm <sup>3</sup>
T1W TSE		
DCE-MRI	↓ Spatial/ ↑ temporal resolution	↑ Spatial/ ↓ temporal resolution
DWI-EPI		↓ Geometric accuracy
Non-EPI DWI	↓↓ SNR	
PseudoCT/MRCAT for MR-only simulation		↓ Geometric accuracy
Treatment simulation, imaging in treatment position	1.5T	3.0T
Brain	Consider for stereotactic treatment	Consider for non-stereotactic treatment
Pelvis: prostate, rectum, cervix, bladder	↓ Resolution	↑ Resolution
Head and neck	↑ Robustness	↑ Resolution
Thorax: esophagus		Artifacts due to motion sensitivity and susceptibility
Abdomen: pancreas, liver		Motion sensitivity concerns
Breast		Susceptibility artifacts
Tumor visualization and response monitoring	1.5T	3.0T
Brain	↓ Contrast, ↓ Resolution	↑ Contrast, ↑ Resolution
Pelvis: prostate, rectum, cervix, bladder	↓ Resolution	
Head and neck		
Thorax: esophagus		Artifacts due to motion sensitivity and susceptibility
Abdomen: pancreas, liver		Motion sensitivity
Breast	DWI better than on 3.0T, ↓ Resolution	↑ Resolution
System-related considerations	1.5T	3.0T
Geometric accuracy		Larger susceptibility changes
SAR		SAR limits in high temporal imaging
Staff	1.5T	3.0T
Expertise	Modest level	High level
Training	Modest level	High level

Not feasible    Not preferable, but possible    Effective, but with some tradeoffs    Most effective



© 2018 Koninklijke Philips N.V. All rights reserved. Specifications are subject to change without notice. Trademarks are the property of Koninklijke Philips N.V. or their respective owners.

4522 991 41071 \* NOV 2018

**How to reach us**

Please visit [www.philips.com/mrirt](http://www.philips.com/mrirt)  
[www.philips.com/radiationoncology](http://www.philips.com/radiationoncology)