



Syntegra

Automated image registration algorithms

Peter Rösch, Ph.D., Thomas Netsch, Ph.D.,
Todd McNutt, Ph.D., James Shoenbill, M.Sc., and
Patricia Root, Product Manager

Rapid technical advances in three-dimensional imaging (e.g., SPECT, MCD, PET, CT, and MRI) have increased the need for relating functional to morphological and anatomical information. Overlaying images in the same spatial orientations facilitates the interpretation of two or more image sets. The aim of image registration is to find a common coordinate system between two image sets in order to correlate a patient's anatomical and functional information.

Syntegra is a multi-modality image registration application for the Philips PET systems and the Philips Pinnacle³® radiotherapy treatment planning system. Syntegra provides manual and point-based image registration, and three automated methods of gray value-based image registration.

Generally, gray value-based methods determine the registration transformation by iteratively optimizing a similarity measure that is calculated from the gray values of both images. The following gray value-based methods offer many advantages when compared to point-based and surface-based methods. Their principal advantages are sub-voxel registration

accuracy and the ability to operate automatically without user interaction.

The three auto-registration methods are defined by the similarity measure that each method uses: cross correlation, local correlation, and mutual information.

Cross Correlation

The Cross Correlation (CC) method seeks to minimize the mean-square gray-value difference between two image sets [1]. In some cases, the images to be registered may not match in gray value but may have a linear relationship between corresponding gray values in the images. Therefore, the CC similarity measure is defined as

$$C^2(M) = \frac{\left(\sum_{i \in M} (I_1(i) - \bar{I}_1)(I_2(i) - \bar{I}_2) \right)^2}{\sum_{i \in M} (I_1(i) - \bar{I}_1)^2 \sum_{i \in M} (I_2(i) - \bar{I}_2)^2}$$

where M is the region of overlap between the two image sets, I_1 and I_2 are the two image sets, i is the index into the image sets, and \bar{I}_1 and \bar{I}_2 are the mean gray value in each image set.

Maximum squared cross correlation corresponds to the optimum alignment only in the presence of a global linear relationship between image gray values of reference and secondary images. For this reason, CC may be utilized when registering two single-modality image sets such as CT/CT or MR/MR using similar protocols.

CC's linear relationship between gray values across the full extent of the regions used in registration is sufficient when both images are from the same image modality and the linear relationship exists across the full range of gray values. When two images are from separate imaging modalities but represent the same information (e.g., MRI and CT both represent anatomy), it is possible to investigate the CC in local sub-regions of the images using local correlation.

Local Correlation

Local Correlation (LC) assumes a local linear relationship between gray values. For example, although the contrast of bone and soft tissue boundaries in MR and CT images is inverted, a local linear relationship (negative slope) exists, at least close to the edges. The assumption is that in small local sub-regions of the image sets, a linear relationship between gray values in the two image sets exists. However, this is dependent on the size and position of the sub-region. If the sub-region is too large, it may intersect with more than two distinct parts of the human anatomy. In this case, the linear relationship approximation may break down.

The LC algorithm applies cross correlation to many small neighborhoods within the image sets instead of to the two image sets overall [2, 3]. The LC similarity measure is thus the sum of the cross correlations of all of the local sub-regions used in the registration, defined as

$$LC^2(M) = \frac{1}{N} \sum_{S_j \in M} C^2(S_j)$$

where C is the cross correlation coefficient for the j th sub-region S_j , and N is the number of sub-regions contained in M .

Local sub-regions are chosen by an evaluation of the voxels within an image set. Sub-regions are determined and chosen as appropriate for the registration process based on the magnitude of variance of image gray values within each sub-region. For this reason, the LC method applies increased importance to high gradient edges present in the image sets.

The LC similarity measure is attractive because it can be formulated as a least-squares problem, which allows for use of faster optimization algorithms. These faster algorithms offer substantial speed improvements (10-20 times faster) over those used in mutual information.

Mutual Information

Mutual Information (MI) [4, 5] may be utilized for a variety of multi-modality automated registration tasks. Probability distributions of the gray values in each image set are calculated and used by the MI similarity equation. The MI algorithm does not rely on an explicit relationship between the gray values in the image sets. Therefore, MI is more suitable in cases of large misalignment or when one of the images to be registered has limited anatomical detail.

MI and normalized MI are based on the gray value histogram, defined as

$$MI = -\sum_{j,k} \left(\frac{P_{j,k}^{2D}}{V} \log \frac{P_{j,k}^{2D}}{P_j^r P_k^f} \right)$$

where V denotes the volume of overlap between the images, and P_j^r and P_k^f are the probabilities of gray value j and k in the reference and secondary image respectively. $P_{j,k}^{2D}$ is the probability that gray values j and k occur in the reference image and at the corresponding position of the secondary image. The probabilities are derived from the bin counts of a joint gray value histogram of the image sets that does not contain spatial information.

Automated registration in Syntegra

The three similarity measures (CC, LC, and MI) provide the capability to automatically register a

broad range of imaging modalities in order to align functional and anatomical information for specific patients. Syntegra provides these methods with appropriate settings for common cross-modality registration procedures (e.g., CT-PET, CT-MRI, etc.) in addition to manual and point-based registration methods.

References

- [1] J.B. Maintz and M.A. Viergever. A survey of medical image registration. *Medical Image Analysis* 2(1): 1-36, 1998.
- [2] P. Rösch, T. Blaffert, and J. Weese. Multi-modality image registration using local correlation. In *Computer Assisted Radiology and Surgery (CARS '99)*, edited by H.U. Lemke, et al. Amsterdam: Elsevier, 1999: 228-232.
- [3] T. Netsch, P. Rösch, A. van Muiswinkel, and J. Weese. Towards real-time multi-modality 3-D medical image registration. *8th International Conference on Computer Vision (ICCV)*. Vancouver, B.C.: IEEE Press, 501-508, 2001.
- [4] F. Maes, A. Collignon, D. Vandermeulen, G. Marchal, and P. Suetens. Multi-modality image registration by maximization of mutual information. *IEEE Transactions on Medical Imaging* 16(2): 187-198, 1997.
- [5] P. Viola and W. M. Wells. Alignment by maximization of mutual information. *International Journal of Computer Vision* 24(2): 137-154, 1997.

Philips Medical Systems is part
of Royal Philips Electronics

INTERESTED?

Would you like to know more about our imaginative products? Please do not hesitate to contact us. We would be happy to provide specific information about our products and services, or put you on our mailing list for news about new product developments, upcoming events or for our clinical journal, Medicamundi. We would be glad to hear from you.

On the web

www.medical.philips.com

Via e-mail

medical@philips.com

By fax

+31 40 27 64 887

By postal service

Philips Medical Systems

Global Information Center

I.B.R.S. / C.C.R.I. Numéro 11088

5600 VC Eindhoven

Pays-Bas / The Netherlands

(no stamp required)

Asia

Regional Office Asia Pacific

Tel: +852 2821 5888

Fax: +852 2527 6727

Europe

Regional Office Netherlands

Tel: +31 40 27 62 614

Fax: +31 40 27 64 250

Germany / Eastern Europe

Regional Office Germany

Tel: +49 1805 767 222

Fax: +49 1805 767 229

Middle East / Africa

Regional Office Dubai

Tel: +971 4 3095 247

Fax: +971 4 3095 200

Latin America

Regional Office Brazil

Tel: +55 11 5188 0764

Fax: +55 11 5188 0761

North America

Regional Office USA

Tel: +1 425 487 7000

Fax: +1 425 487 8130

© Royal Philips Electronics N.V. 2003

All rights are reserved. Reproduction in whole or in part is prohibited without the prior written consent of the copyright holder.

Philips Medical Systems Nederland B.V. reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

Printed in the USA.
4535 983 02478/870 10M 04/03

