



Product name	XW6600 - 8999 Computer (Reconstructor / Ambient Experience / Host)
Identification code(s)	XW6600: 4522 132 7684x, 4522 132 7685x, 4522 132 7686x , 4522 133 0044x XW8600: 4522 132 7694x, 4522 132 7695x
Total weight [kg]	15 kg – 25 kg
Manufacturer	Hewlett Packard
Address	Orteliuslaan 1000
Zip code	3528 BD Utrecht
Country	The Netherlands
Electronic info	452213276840_00_5790_Computer_XW6600-8999_v02.doc

Recycling information

Batteries	Items, type	Location
	Lithium battery type CR2032 (3V)	Figure 1a

Hazardous substances	Items	Location
	Printed Circuit Boards maximal 8 add-in cards (memory) with a surface greater than 10 square cm.	Figure 1b
	Printed Circuit Boards system board (motherboard) with a surface greater than 10 square cm.	Figure 1b
	Printed Circuit Boards maximal 7 add-in cards (PCI) with a surface greater than 10 square cm	Figure 2
	Printed Circuit Boards two in the PSU with a surface greater than 10 square cm.	Figure 3a and 3b
	Four attached to the power supply Printed Circuit Board cards. Electrolytic capacitors / condensers that measure greater than 2.5 cm in diameter or height.	Figure 3a and 3b

Name: Recycling Passport Computers XW6600-8999		12NC: 4522 132 76840		Version: 00 Doc.Date: 2008-04-23	
Author: N.K. Tran		Pag.: 1 of 4	Group No: 5790	Section: -	Status: Authorized

Locations as mentioned in the Passport (Information in pictures)




Figure 1a removing the system battery



Press on the release tab of the battery holder. Rotate the battery beyond the latch and lift it straight up.




Figure 1b removing the system board and DIMMs



Slide the board forward (1) to disengage the metal mounting standoffs from the chassis. Lift the board out of the chassis (2).

Name: Recycling Passport Computers XW6600-8999	12NC: 4522 132 76840	Version: 00 Doc.Date: 2008-04-23
Author: N.K. Tran	Pag.: 2 of 4	Group No: 5790 Section: - Status: Authorized




Figure 2 removing the graphics- and any other expansion cards

Open the PCI card retention clamp. Press the retention clamp levers down to unsnap, and then rotate them up (1). Pull the PCI card up and out of the chassis (2).




Figure 3a removing the power supply from the chassis

Disconnect or cut electrical connections to the power supply. Remove the four screws securing the power supply to the workstation chassis (1). Slide the power supply toward the front of the chassis and lift it out of the chassis (3).

Name: Recycling Passport Computers XW6600-8999		12NC: 4522 132 76840		Version: 00 Doc.Date: 2008-04-23	
Author: N.K. Tran		Pag.: 3 of 4	Group No: 5790	Section: -	Status: Authorized




Figure 3b removing the desired components from both Printed Circuit Board cards

Ensure that all capacitors are safely discharged. Use a diagonal cutter to remove the desired components from both Printed Circuit Board cards. (Note that some components may need to be unsoldered if a diagonal cutter does not fit between the component and the Printed Circuit Board Card.)

Recycling Passport Revision History

Revision	Description	Date
00	Initial version	2008-03-18
00	Merge of recycling passports XW6600, XW8000, XW8200	2008-04-22
00	Updated after review HW team	2008-04-23

Recycling Passport Approval Signatures

Written By N.K. Tran	Reviewed By Geert-Jan Plattèl
Approved By N.K. Tran	Approved By Geert-Jan Plattèl

Name: Recycling Passport Computers XW6600-8999		12NC: 4522 132 76840		Version: 00 Doc.Date: 2008-04-23	
Author: N.K. Tran		Pag.: 4 of 4	Group No: 5790	Section: -	Status: Authorized