





PHILIPS HEALTHCARE	XL+ Product Recycling Passport	Page 1 of 6
-------------------------------	---------------------------------------	-------------

Product name:	Heartstart XL+	
Identification code(s)	861290	
Total weight (in kg)	Up to 9 kg	
Producer/ Manufacturer	Name company:	Philips Healthcare
	Address:	3000 Minuteman Road, Andover, MA
	Zip code:	01810
	Country:	USA

Recycle Info	Items:	Location
<p>Special attention</p> 	<p>WARNING: Remove all power sources (unplug AC power and remove the external battery, see Figure 1). Make sure the defibrillator is disarmed before disassembling the defibrillator otherwise dangerous voltages may be present on components and connections. After opening the defibrillator make sure the high voltage capacitor is discharged by using the Defibrillator Discharge Tool, Philips part number M2475-69573. Connect the hooked end of the discharge tool to wire loop #1 on the Therapy PCA, see figure 2. Then touch the pointed end of the discharge tool to wire loop #2. Hold the discharge tool in place for at least 10 seconds. Failure to discharge the defibrillator could result in serious injury or death.</p> <p>To order the Defibrillator Discharge Tool contact your local Philips Healthcare Sales Office or your authorized Philips Healthcare Dealer or Distributor</p>	Figures 1 & 2

	Items:	Location
	None	

<p>PHILIPS HEALTHCARE</p>	<p align="center">XL+ Product Recycling Passport</p>	<p align="right">Page 2 of 6</p>
--------------------------------------	---	----------------------------------

<p>Batteries</p>  <p>To be Removed</p>	<table border="1"> <thead> <tr> <th align="center" colspan="2">Type:</th> <th align="center">Location</th> </tr> </thead> <tbody> <tr> <td>External Lithium Ion</td> <td></td> <td>Figure 1</td> </tr> <tr> <td>Internal Lithium</td> <td></td> <td>Figure 3</td> </tr> <tr> <td></td> <td></td> <td></td> </tr> <tr> <td></td> <td></td> <td></td> </tr> <tr> <td></td> <td></td> <td></td> </tr> <tr> <td></td> <td></td> <td></td> </tr> </tbody> </table>	Type:		Location	External Lithium Ion		Figure 1	Internal Lithium		Figure 3																
Type:		Location																								
External Lithium Ion		Figure 1																								
Internal Lithium		Figure 3																								
<p>Hazardous</p>  <p>To be Removed</p>	<table border="1"> <thead> <tr> <th align="center" colspan="2">Substances:</th> <th align="center">Location</th> </tr> </thead> <tbody> <tr> <td>LCD display backlight tube (CCFL) contains mercury (Hg)</td> <td></td> <td>Figure 4</td> </tr> <tr> <td>Many of the printed circuit boards and wire solder connections contain leaded solder</td> <td></td> <td>none</td> </tr> <tr> <td></td> <td></td> <td></td> </tr> <tr> <td></td> <td></td> <td></td> </tr> <tr> <td></td> <td></td> <td></td> </tr> <tr> <td></td> <td></td> <td></td> </tr> <tr> <td></td> <td></td> <td></td> </tr> </tbody> </table>	Substances:		Location	LCD display backlight tube (CCFL) contains mercury (Hg)		Figure 4	Many of the printed circuit boards and wire solder connections contain leaded solder		none																
Substances:		Location																								
LCD display backlight tube (CCFL) contains mercury (Hg)		Figure 4																								
Many of the printed circuit boards and wire solder connections contain leaded solder		none																								




Figure 1




Figure 2




Figure 3

PHILIPS HEALTHCARE	XL+ Product Recycling Passport	Page 6 of 6
-------------------------------	---------------------------------------	-------------




Figure 4