



# Bispectral Index (BIS)

Monitoring level of consciousness with the compact BISx Power Link™

Philips Bispectral Index (BIS®) measurement uses Aspect Medical Systems' XP platform technology. It provides depth of consciousness and sedation monitoring for use in the OR, ICU, and other clinical settings. A sensor on the patients forehead translates brain electrical activity into a single BIS value, ranging from 0 to 100, with 100 signifying full awareness and 0 indicating no brain activity.

The Philips BIS module provides numerics, graphic trends, and high-resolution trends for the following values:

- Bispectral Index®
- Electromyographic strength (EMG)
- Signal quality indicator
- Suppression ratio
- EEG waveform
- Total power
- Spectral Edge Frequency (SEF)
- Burst (available only via Extend sensors)

BIS measurements are resistant to artifact from electrocautery devices. The algorithm detects and filters interference from EMG.

## Benefits for anesthesia and critical care


BIS measurements can help clinicians formulate the precise type and optimal dosages of anesthetic or sedative medication for each patient.

BIS monitoring is most widely used in the OR to help:

- Regulate anesthetic drug use
- Decrease the incidence of post-operative side effects such as nausea and vomiting
- Reduce length of stay in the PACU (recovery room)
- Prevent intraoperative awareness

BIS monitoring can also be useful during outpatient surgery and conscious sedation.

BIS measurements can be arranged on the screen to suit patient requirements and user preferences.



**PHILIPS**

## Compatibility

BIS information is displayed alongside other key parameters on these Philips patient monitors:

- IntelliVue MP20 and higher with IntelliVue release C and higher
- CMS 2002
- V24/V26

Please ask your sales representative for details on compatibility.

In the ICU, BIS monitoring has been shown to reduce recall of unpleasant experiences and provide objective sedation assessment during:

- Mechanical ventilation
- Neuromuscular blockade
- Barbiturate coma
- Bedside procedures

Well managed sedation levels in the ICU can also aid in ventilator weaning.

## Greater awareness of patient-specific reactions to anesthesia and sedation

Patients respond differently to anesthesia, and BIS monitoring can be particularly valuable in patients whose responses are least predictable, including:

- Geriatric and pediatric patients
- Frail patients
- Patients undergoing lengthy surgeries
- Patients with conditions known to affect responsiveness to anesthesia

## BISx Power Link works with a range of sensors

The plug-and-play BISx Power Link™ is compact and portable. It has an integrated clamp for quick and easy mounting and connects seamlessly by cable to the Philips BIS module. No extra mounts are needed. BIS sensors are tailored to the special requirements of ICU, OR, and pediatric patients.



BIS is one of many Philips measurements and technologies designed for the demands of anesthesia care.



The BIS module and compact, portable BISx Power Link for software and signal processing.

M1034Ax	BISx Power Link™
M1034AX #A01	BISx Power Link™ plus BIS Module
M1034AX #A02	BISx Power Link™ plus BIS interface
M1034A	BIS module only
M8001A #J16	BIS interface board for MP20
M8002A #J16	BIS interface board for MP30

## BIS sensors:

Quatro Sensor has four electrodes that provide full BIS XP performance in the OR, ICU, and other clinical settings

Extend Sensor has four electrodes that provide full BIS XP performance, especially designed for the ICU (currently available in the US only)

Pediatric Sensor has four electrodes suited to the small size of pediatric patients\*

Semi Reusable Sensor (SRS) Kit\*\* consists of 100 disposable sensors and a cable that is reusable up to 100 times\*

\*\* Not available in the US or Japan.

+ Supported only on patient monitors with release C and higher.

## Philips Commitment to Measurement Technologies

Philips is committed to providing best-in-class standard clinical measurements as well as innovative measurements to support clinicians' decisions at the patient's side.

- Maintaining and advancing the performance of existing, widely used standard-of-care measurements
- Investing heavily in research, development, and clinical validation of new, innovative parameters and algorithms
- Working with strategic partners to integrate next-generation measurements and technologies
- Providing interfaces to more than 100 third-party specialty measurement devices through the Philips VueLink module



© 2006 Koninklijke Philips Electronics N.V.  
All rights are reserved.

Philips Medical Systems Nederland B.V. reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

Bispectral Index, BIS, BISx Power Link, and Aspect are trademarks or registered trademarks of Aspect Medical Systems.

Philips Medical Systems is part of Royal Philips Electronics

www.medical.philips.com  
medical@philips.com  
fax: +31 40 27 64 887

Printed in The Netherlands  
4522 962 13641/862 \* APR 2006

Philips Medical Systems  
Global Information Center  
P.O. Box 1286  
5602 BG Eindhoven  
The Netherlands