



**DICOM 3.0
Conformance Statement**

For

**(WPS)
Witt PACS Server v1.0**

**Revision 1.5
August 11, 2004**

*© Witt Biomedical 2001
All rights reserved*

TABLE OF CONTENTS

1. INTRODUCTION	6
1.1 Scope.....	6
1.2 Content Structure	6
1.3 Intended Audience	6
1.4 Requirements and Use.....	6
1.5 Acronyms and abbreviations	7
1.6 References.....	7
2. IMPLEMENTATION MODEL VERIFICATION	8
3. IMPLEMENTATION MODEL.....	8
3.1 Application Data Flow Diagram.....	8
3.1.1 Verification Application Entity.....	8
3.1.1.1 Presentation Context Table	8
3.1.2 Image Storage Application Entity.....	9
3.1.2.1 Presentation Context Table	9
3.1.3.....	11
3.1.4 Query/Retrieve Application Entity	12
Accepted Q/R SOP Class.....	12
3.2 Functional Definition of AE	13
3.2.1 Verification Application Entity.....	13
3.2.2 Image Storage Application Entity.....	13
3.2.3 Q/R Application Entity	14
4. COMMUNICATION PROFILES.....	15
4.1 Supported Protocol Stacks (parts 8, 9)	15
4.2 TCP/IP Stack.....	15
4.3 Physical Media Support	15
5. CONFIGURATION	15
5.1 Verification AE.....	15

5.2	Image Storage AE	15
5.3	Query/Retrieve AE.....	16
6.	SUPPORT OF BASIC/EXTENDED CHARACTER SETS.....	16

Revision History

Revision	Date	Author	Reason for Change
1.0	8/17/01	CM	Preliminary version
1.1	12/17/01	CM	Module addition
1.2	1/21/01	CM	Presentation Syntax Table
1.3	6/23/03	CM	Typographical adjustment
1.4	10/15/03	CM	Addition to Accepted Q/R SOP Class
1.5	8/11/04	JC	Addition to Accepted Storage SOP Class and Q/R SOP Class

Copyright Statement

Witt Biomedical has used every resource available to ensure the accuracy of this document. Although we have taken all those precautions Witt Biomedical assumes no liability for errors or omissions and reserves the right to make changes without further notice to any products herein to improve reliability, function, or design.

Witt Biomedical Corporation provides this guide without warranty of any kind, implied or expressed, including, but not limited to, the implied warranties of merchantability and fitness for a particular purpose.

Witt Biomedical Corporation may make improvements and/or changes in the product(s) and/or module(s) described in this document at anytime. This document contains proprietary information, which is protected by copyright (All rights are reserved).

This document may not be reproduced, photocopied, converted or translated to another language without prior written permission from Witt Biomedical Corporation. Witt Biomedical has to right to make changes to this document or revise it without any obligation to provide notification of such revision or changes.

Disclaimer

Witt Biomedical certifies that the WPS software is in compliance with the ACR-NEMA DICOM 3.0 standard. However, due to the flexibility of DICOM, the user and/or company must perform integration testing to verify that the WPS is compatible and meets the requirements of the integration with another system. The testing must include but not be limited to sending DICOM Datasets (images/video) for all type of images you wish to

transfer. If the results of your testing are not conclusive please contact your representative to guide you through your testing process.

1. Introduction

1.1 Scope

This document states the conformance statement of Witt Biomedical DICOM compliant software/hardware family to DICOM 3.0 standard. It applies to:

Witt PACS Server (Hereinafter WPS):

- Verification & Store Module (Included)
- Query & Retrieve Module (Optional)

Please note that each module has its own capabilities, which are identified when considered appropriate.

1.2 Content Structure

The DICOM conformance statement consists of sections 2 through 7. It follows the content requirements of DICOM PS 3.2

1.3 Intended Audience

This Conformance Statement is intended for software engineers, system integrators, field engineers, and biomedical technicians. The audience is assumed to have a practical and working knowledge of DICOM standard and software interfaces in general.

1.4 Requirements and Use

Since the DICOM interface option and/or modules are not available if not purchased, System integrators who wish to implement a DICOM compliant hardware device and/or software from another manufacturer, need to contact the appropriate authority to activate or purchase the appropriate options/modules. In some circumstance hardware upgrade might be required to accommodate such options.

1.5 Acronyms and abbreviations

The following acronyms and abbreviations are used in this Conformance Statement:

AE	Application Entity
DICOM	Digital Imaging and Communication in Medicine
FIFO	First In First Out
HL7	Health Level Seven
IP	Internet Protocol
JPEG	Joint Photographic Experts Group (compression format)
SCP	Service Class Provider
SCU	Service Class User
SOP	Service Object Pair
TCP/IP	Transmission Control Protocol/Internet Protocol
UID	Unique Identifier
Q/R	Query and Retrieve
WPS	Witt PACS Server

1.6 References

DICOM Digital Imaging and Communication in Medicine, NEMA PS 3.1-3.14, 1999

2. Implementation Model Verification


In order for foreign applications (SCU) to verify if they can request a connection and establish one they must perform a verification request. WPS Image Storage service is responsible to respond to such request. If such request fails make sure that the foreign AE is configured on the WPS Image Store configuration.

3. Implementation Model

Each DICOM SCU and SCP is an Application Entity (AE).

3.1 Application Data Flow Diagram

3.1.1 Verification Application Entity




3.1.1.1 Presentation Context Table

Any of the presentation context found in the following reference table are accepted for the DICOM Store SCP to receive the verification request (C-ECHO).

Presentation Context Table					
Abstract Syntax		Transfer Syntax		Role	Extended Negotiation
Name	UID	Name	UID		
Verification	1.2.840.10008.1.1	Implicit VR Little Endian	1.2.840.10008.1.2	SCP	None
Verification	1.2.840.10008.1.1	Explicit VR Little Endian	1.2.840.10008.1.2.1	SCP	None

3.1.2 Image Storage Application Entity



3.1.2.1 Presentation Context Table

Here is a list of image storage classes accepted by the Image Storage AE. Any of the presentation contexts in the following table are accepted for the DICOM Store SCP to receive images. The user can restrict accepted presentation context in the Store SCP configuration panel.

Presentation Context					
Abstract Syntax		Transfer Syntax		Role	Extended Negotiation
Name	UID	Name	UID		
Computed Radiography Image Storage	1.2.840.10008.5.1.4.1.1.1	See Table below	See Table below	SCP	None
CT Storage	1.2.840.10008.5.1.4.1.1.2				
MR Image Storage	1.2.840.10008.5.1.4.1.1.4	See Table below	See Table below	SCP	None
Secondary Capture Image Storage	1.2.840.10008.5.1.4.1.1.7	See Table below	See Table below	SCP	None
Ultrasound Multi-frame Image Storage (Retired)	1.2.840.10008.5.1.4.1.1.3	See Table below	See Table below	SCP	None
Ultrasound Multi-frame Image Storage	1.2.840.10008.5.1.4.1.1.3.1	See Table below	See Table below	SCP	None

Ultrasound Image Storage (Retired)	1.2.840.10008.5.1.4.1.1.6	See Table below	See Table below	SCP	None
Ultrasound Image Storage	1.2.840.10008.5.1.4.1.1.6.1	See Table below	See Table below	SCP	None
X-ray Angiographic Bi-plane Image Storage	1.2.840.10008.5.1.4.1.1.12.3	See Table below	See Table below	SCP	None
X-ray Angiographic Image Storage	1.2.840.10008.5.1.4.1.1.12.1	See Table below	See Table below	SCP	None
X-ray RadioFluoroscopic Image Storage	1.2.840.10008.5.1.4.1.1.12.2	See Table below	See Table below	SCP	None
Nuclear Medicine Image Storage	1.2.840.10008.5.1.4.1.1.20	See Table below	See Table below	SCP	None
Nuclear Medicine Image Storage (Retired)	1.2.840.10008.5.1.4.1.1.5	See Table below	See Table below	SCP	None
PET Image Storage (DICOM 98)	1.2.840.10008.5.1.4.1.1.128	See Table below	See Table below	SCP	None
Digital Xray Pres Image Storage	1.2.840.10008.5.1.4.1.1.1.1	See Table below	See Table below	SCP	None
Digital Mammo Pres Image Storage	1.2.840.10008.5.1.4.1.1.1.2	See Table below	See Table below	SCP	None
Digital Intra Oral Pres Image Storage	1.2.840.10008.5.1.4.1.1.1.3	See Table below	See Table below	SCP	None
Visible Light Image Storage	1.2.840.10008.5.1.4.1.1.77.1.1	See Table below	See Table below	SCP	None
Visible Light Micro-Film Image Storage	1.2.840.10008.5.1.4.1.1.77.1.2	See Table below	See Table below	SCP	None
Visible Light Slide Image Storage	1.2.840.10008.5.1.4.1.1.77.1.3	See Table below	See Table below	SCP	None
Visible Light Photo Image Storage	1.2.840.10008.5.1.4.1.1.77.1.4	See Table below	See Table below	SCP	None
Structured Reporting (SR) Detailed Storage	1.2.840.10008.5.1.4.1.1.88.3	See Table below	See Table below	SCP	None
Basic Text SR	1.2.840.10008.5.1.4.1.1.88.11	See Table below	See Table below	SCP	None
Enhanced SR	1.2.840.10008.5.1.4.1.1.88.22	See Table below	See Table below	SCP	None
Comprehensive SR	1.2.840.10008.5.1.4.1.1.88.33	See Table below	See Table below	SCP	None

Each of the storage SOP classes listed above can be transferred using one of the following transfer syntaxes.


Note: If client initiate a multi-file send during a single association, all file must be of the same SOP & Transfer syntax.

Transfer Syntax Name	UID
Implicit VR Little Endian	1.2.840.10008.1.2
Explicit VR Little Endian	1.2.840.10008.1.2.1
Default Lossy JPEG Compressed	1.2.840.10008.1.2.4.50
Default Lossless JPEG Compressed	1.2.840.10008.1.2.4.70
RLE Compressed	1.2.840.10008.1.2.5
Others (Optional)	
Explicit VR Big Endian	1.2.840.10008.1.2.2
JPEG Extended (Process 2 & 4)	1.2.840.10008.1.2.4.51
JPEG Extended (Process 3 & 5)	1.2.840.10008.1.2.4.52
JPEG Spectral Selection, Non-Hierarchical (Process 6 & 8)	1.2.840.10008.1.2.4.53

JPEG Spectral Selection, Non-Hierarchical (Process 7 & 9)	1.2.840.10008.1.2.4.54
JPEG Full Progression, Non-Hierarchical (Process 10 & 12)	1.2.840.10008.1.2.4.55
JPEG Full Progression, Non-Hierarchical (Process 11 & 13)	1.2.840.10008.1.2.4.56
JPEG Lossless, Non-Hierarchical (Process 14)	1.2.840.10008.1.2.4.57
JPEG Lossless, Non-Hierarchical (Process 15)	1.2.840.10008.1.2.4.58
JPEG Extended, Hierarchical (Process 16 & 18)	1.2.840.10008.1.2.4.59
JPEG Extended, Hierarchical (Process 17 & 19)	1.2.840.10008.1.2.4.60
JPEG Spectral Selection, Hierarchical (Process 20 & 22)	1.2.840.10008.1.2.4.61
JPEG Spectral Selection, Hierarchical (Process 21 & 23)	1.2.840.10008.1.2.4.62
JPEG Full Progression, Hierarchical (Process 24 & 26)	1.2.840.10008.1.2.4.63
JPEG Full Progression, Hierarchical (Process 25 & 27)	1.2.840.10008.1.2.4.64
JPEG Lossless, Hierarchical (Process 28)	1.2.840.10008.1.2.4.65
JPEG Lossless, Hierarchical (Process 29)	1.2.840.10008.1.2.4.66

3.1.3

3.1.4 Query/Retrieve Application Entity



Accepted Q/R SOP Class

WPS provides Standard Conformance to the following DICOM V3.0 **Query/Retrieve** SOP Class as an SCP (Server).

Description	SOP Class UID
Patient Root Query/Retrieve C-FIND Request	1.2.840.10008.5.1.4.1.2.1.1
Patient Root Query/Retrieve C-MOVE Request	1.2.840.10008.5.1.4.1.2.1.2
Study Root Query/Retrieve C-FIND Request	1.2.840.10008.5.1.4.1.2.2.1
Study Root Query/Retrieve C-MOVE Request	1.2.840.10008.5.1.4.1.2.2.2
Modality Worklist C-FIND Request	1.2.840.10008.5.1.4.31

3.2 Functional Definition of AE

3.2.1 Verification Application Entity

The Verification AE is a SCP module that permanently waits for a DICOM 3.0 verification request from a remote Verification AE. Upon connection it returns a response to the remote Verification AE based on the WPS current configuration.

The response can be one of the following:

- Success (Remote AE Accepted)
- Rejected (Remote AE Rejected)

The response is based on an internal query of acceptable remote AE connections.

3.2.2 Image Storage Application Entity

The Image Storage AE uses DICOM 3.0 Storage service (SCP) to store single and multi-frame image/video. Supported image/video SOP CLASS are described in section 3.1.2.

WPS Image Storage service can accept a limited amount of simultaneous connections. The number of simultaneous connections is defined in the main Image Storage SCP configuration screen.

WPS can accept all AE Titles SCU (in test mode) or accept a limited list of AE for security purposes.

WPS will accept and store the image “as-is” without modification of its type (RLE, JPEG...) or attributes.

The user can choose which transfer syntax they wish to accept for each SOP Class. This is done in the Image Storage configuration screen.

Upon reception of a DICOM file (single or multi-frame), WPS will verify if the image was already received previously and process accordingly.

Each DICOM file will be placed on the WITT Network and kept on the server (including long term storage) ready to be used for review by WITT review stations or be queried by foreign review stations.

3.2.3 Q/R Application Entity

The Query and Retrieve AE uses DICOM 3.0 Q/R service (SCP) to query stored studies and DICOM 3.0 Store (SCU) to return selected images. Supported image/video SOP CLASS are described in section 3.1.3.

WPS Q/R Service can accept a limited quantity of simultaneous connections. The number of simultaneous connections is defined in the Q/R main configuration page.

When WPS Q/R service receives a C-FIND command, it initiates a query on its own or remote SQL study server. WPS will base the query on the type of C-FIND command received. Once the local query is executed the information is returned in DICOM format to the Q/R SCU.

From that point the remote Q/R SCU can refine its search or send a C-MOVE command to the Q/R SCP. WPS Q/R SCP receives the appropriate command and appends a record to the STORE SCU queue, which will send the requested file to the specified STORE SCP.

4. Communication Profiles

4.1 Supported Protocol Stacks (parts 8, 9)

Witt PACS network apparatus are all using DICOM upper layer protocol as defined in Parts 8 and 9 of DICOM standard.

Our system is using TCP/IP stack on all DICOM compliant devices including print. (Excluding physical media transfer on CD-R, DVD, MO)

4.2 TCP/IP Stack

Witt PACS devices are all using TCP/IP stack via Microsoft Windows Winsock interface.

4.3 Physical Media Support

Witt Biomedical recommends using at least 100BASE-T (IEEE 802.3) network and network devices. Slower network, such as 10BASE T, would provide unacceptable user response times for almost all modality data sets.

5. Configuration

Every AE can configure the following parameters:

- Acceptable foreign SCU (IP Address and/or AE Title)
- Port number for SCU & SCP
- Number of simultaneous connection per service
- Proprietary Network Information

5.1 Verification AE

5.2 Image Storage AE

- Network directory structure (if storing on foreign server... (NFS, FTP of Microsoft Network))
- Location and type of SQL Server in use (Microsoft, Oracle, MySQL)
- Acceptable SOP CLASS
- Acceptable Transfer Syntax for each SOP CLASS
- Long Term Storage Configuration

5.3 Query/Retrieve AE

- Acceptable Storage SCP Redirection Port and AE Title

6. Support of Basic/Extended Character Sets

WPS supports the following character set: ISO-IR (100) Latin alphabet #1