

FieldStrength

Publication for the Philips MRI Community

Issue 41 – September 2010



mDIXON being developed to simplify and accelerate liver MRI

The Philips mDIXON scan is being designed to create four image types in one breath hold



This article is part of
Field Strength issue 41,
September 2010

PHILIPS
sense and simplicity

mDIXON being developed to simplify and accelerate liver MRI

The Philips mDIXON scan is being designed to create four image types in one breath hold


A liver MRI exam typically includes T2-weighted imaging, chemical shift imaging, pre-contrast T1-weighted, and dynamic post-contrast imaging. Currently, three separate scans are required to create these four image datasets. The mDIXON sequence is being developed to provide four image types in one breath hold: water and fat images as well as in-phase and opposed-phase images.


The mDIXON sequence is being developed to provide the same images as the pre-contrast and chemical shift sequences, in only one single scan that fits in a breath hold of about 20 seconds. Because this mDIXON scan is designed to be faster, it allows using higher in-plane and through-plane spatial resolution within one breath hold. Shahid Hussain, MD, PhD, University of Nebraska Medical Center, USA and Thomas Perkins, PhD, Senior Clinical Scientist at Philips, are cooperating to optimize the mDIXON sequence for liver MRI.

Philips mDIXON differs from previous Dixon techniques

The original Dixon technique requires a time-consuming acquisition of in-phase and opposed-phase gradient echo images. The in-phase and opposed-phase images are then added together to get water-only images, and subtracted to get fat-only images. These imaging sequences have limited possibilities to simultaneously optimize for spatial resolution, slice thickness, and scan time, and still fit in a single breath hold in the case of liver imaging. The Philips mDIXON sequence is based on a unique reconstruction algorithm, offering

Schematic comparison of mDIXON and standard methods.





freedom in echo time setting, so that it does allow for optimization of key imaging parameters and still fits in one breath hold. Arbitrary echo times are used to generate reconstructed water-only images, fat-only images, in-phase images and opposed-phase images.

Time-consuming fat-suppression RF pulses are not needed, and this gain in speed allows to select better spatial resolution for mDIXON than for e-THRIVE within the same breath hold time. In addition, mDIXON performs a special optimization designed to obtain more uniform fat suppression than with e-THRIVE. On 3.0T systems, a two-point mDIXON sequence may be ideal, because at higher field strengths, the frequency difference between fat and water increases, shortening the fat-water amplitude modulation period such that the two acquired mDIXON echo times can be short, thus enabling a higher data sampling rate.

References

1. S Hussain, T Perkins
Preliminary Clinical Experience with a Multiecho 2-Point DIXON (mDIXON) Sequence at 3T as an Efficient Alternative for Both the SAR-Intensive Acquired In- And Out-Of-Phase Chemical Shift Imaging as Well as for 3D Fat-Suppressed T1-Weighted Sequences Used for Dynamic Gadolinium-enhanced Imaging
Proc. Intl. Soc. Mag. Reson. Med, 18:556 (2010)
2. WT Dixon
Simple Proton Spectroscopic Imaging
Radiology, 153:189-194 (1984)
3. H Eggers, B Brendel, G Herigault
Comparison of Dixon Methods for Fat Suppression in Single Breath-Hold 3D Gradient-Echo Abdominal MRI
Proc. Intl. Soc. Mag. Reson. Med., 17:2705 (2009)
4. S Hussain
Liver MRI: Correlation with Other Modalities and Histopathology
Springer-Verlag, Berlin Heidelberg 2007.