

# Designed with you and your patients in mind

The OmniLab Advanced + multi-mode titration system





With the most modes of any titration device, a quiet platform, customizable reports, user-friendly features, and value-added benefits, you'll find the OmniLab Advanced + to be an ideal system for satisfying your lab titration needs.

# Your ally in the lab

We back our commitment to help you achieve lab efficiency and successful titration studies by offering the OmniLab Advanced + titration system. Combined with our OmniLab Direct and OmniLab Connect, the OmniLab Advanced + offers enhanced features designed to help your staff perform titration studies, from beginning to end, easily and efficiently, on even your most complicated patients.

#### Nine therapy modes for enhanced flexibility

| CPAP modes  |  |
|---|--|
| CPAP  | For OSA patients   |
| Auto CPAP   | For OSA patients   |
| Bi-level modes  |  |
| Auto Bi-level (with Bi-Flex)  | For OSA and CPAP rescue patients   |
| S (with AVAPS)  | For OSA patients   |
| S/T (with AVAPS)  | For COPD, obesity hypoventilation, and neuromuscular disorder patients   |
| autoSV (with Bi-Flex)   | For Complex patients with central or mixed apneas and periodic breathing |
| T (with AVAPS)  | For COPD and neuromuscular disorder patients                             |
| PC (with AVAPS)   | For neuromuscular disorder patients (ALS)                                |
| AVAPS-AE mode   |  |
| AVAPS-AE (found in Trilogy series ventilators and select International devices) | For respiratory insufficiency patients (COPD, OHS and NMD)               |

Along with our Advanced therapy options—autoSV and AVAPS—clinicians and sleep technicians can now use AVAPS-AE to assist them with titrating their most complicated patients.

AVAPS-AE, supported by Trilogy ventilators, allows:

- AVAPS to maintain a targeted tidal volume
- Auto EPAP to maintain patent upper airway at a comfortable pressure
- Auto backup rate to target the patient's spontaneous respiratory rate
- Adjustable AVAPS to allow the clinician to control the speed at which pressure will change

# Optimize your sleep lab workflow

From the market leaders in lab titration systems, we now bring you the OmniLab Advanced + system, with OmniLab Direct and OmniLab Connect. Together they provide you with a cohesive, harmonized solution to help you and your sleep technicians operate efficiently, easily and intuitively. The simplicity and ease of use of the titration system, including nine therapy modes, give you more options and flexibility to titrate more types of patients than ever before.

### OmniLab Advanced + advantages

Nine therapy modes—the most available in the industry—to provide you with range and flexibility for all of your titration studies, including pediatric patients.\*

#### Easy-to-navigate color screen—

for use at bedside when setting up your patients to facilitate settings and patient comfort.

Minimal noise output—quiet operation.



#### Detailed titration settings report-

time stamps patient settings and pressure changes directly into a summary report that you can edit and print for later follow up and physician review.

Prescription report—captures the patient's final therapeutic pressure settings into a customizable, printed prescription to help save additional processing time.

Heated Tube humidifier option—for enhanced humidification and reduction of rainout to facilitate patient comfort and compliance.

OmniLab Direct software—allows you to intuitively and remotely control the OmniLab Advanced + device from your control room for easy therapy changes and minimal patient disruptions.

#### OmniLab Connect—provides

flexibility to use up to eight channels when interfacing between OmniLab Advanced + and a non-Philips Respironics' PSG system. OmniLab Connect features an updated setup screen that allows you to create mode-specific montages on your PSG system and select from 24 signals, for specific modes, such as autoSV.

Mask leak detection—displays pertinent mask information, including total and unintentional leak. This information is displayed prior to the beginning of titration, and throughout the night, to help technicians quickly troubleshoot any possible mask issues.

Access to My.Respironics.com—for convenient access to software upgrades and enhancements.

Software upgrades and product notices available—through our SleepDx registration website.


<sup>\*</sup> The effectiveness of the A-Flex and C-Flex+ therapy features, and the autoSV mode has not been established for pediatric patients at this time.



## Settings and reports at your fingertips



OmniLab Direct allows patient settings and modes to be changed remotely from the control room



Prescription tool that captures final therapeutic settings into an automatically generated and customizable prescription

# Specifications and ordering information

| 8 75" (I) × 7 25" (w) × 4 25" (b)   |   |
|---|---|
| , ,   |   |
| · · · · · · · · · · · · · · · · · · ·   |   |
| ,   |   |
| ·   |   |
| •   |   |
|   |   |
| Integrated heated humidification, Heated Tube humidification,                           |   |
|   |   |
| lighted color screen, automatic on/off, mask off alert                                  |   |
| 4 20 11 0   |   |
| 2   |   |
| 4   |   |
| 4   | AP 4 – 30 cm H <sub>2</sub> O (except with Bi-Flex turned on,   |
| ,   | 10.4.05   |
| EPAP 4 – 25 cm H <sub>2</sub> O, IPAP 4 – 25 cm H <sub>2</sub> O                        |   |
| EPAP 4 – 25 cm H <sub>2</sub> O, IPAP 4 – 30 cm H <sub>2</sub> O                        |   |
| EPAP $4 - 30 \text{ cm H}_2\text{O}$ , Pressure Support $0 - 26 \text{ cm H}_2\text{O}$ |   |
| $4-30 \text{ cm H}_2\text{O}$ , Pressure Support $2-26 \text{ cm H}_2\text{O}$          |   |
|   |   |
|   | Part number   |
| ted Tube humidifier   | 1111124   |
|   | 1111551   |
| 50 and 60 Series lids)  | 1063785   |
|   | SYS1HT15  |
|   | 1111076   |
|   | Log in to My.Respironics.com or   |
|   | use PN 1112033 to receive a CD  |
|   | 4 – 20 cm H <sub>2</sub> O<br>4 – 20 cm H <sub>2</sub> O<br>4 – 20 cm H <sub>2</sub> O<br>EPAP 4 – 25 cm H <sub>2</sub> O, IPA<br>then 4-25)<br>EPAP 4 – 25 cm H <sub>2</sub> O, IPA<br>EPAP 4 – 25 cm H <sub>2</sub> O, IPA<br>EPAP 4 – 30 cm H <sub>2</sub> O, Pro- |



We cover your lab diagnostic equipment and titration needs





Philips Respironics is working in harmony with the care team to produce solutions that can lead to healthier patients, healthier practices and healthier businesses.

### Philips Healthcare is part of Royal Philips

How to reach us

www.philips.com/healthcare healthcare@philips.com

Asia

+49 7031 463 2254

Europe, Middle East, Africa +49 7031 463 2254

Latin America +55 11 2125 0744

North America +1 425 487 7000 800 285 5585 (toll free, US only) Philips Respironics 1010 Murry Ridge Lane Murrysville, PA 15668

Customer Service +1 724 387 4000

800 345 6443 (toll free, US only)

Philips Respironics Asia Pacific

+65 6882 5282

Philips Respironics Australia +61 (2) 9947 0440

1300 766 488 (toll free Australia only)

Philips Respironics China +86 400 828 6665 +86 800 828 6665 Philips Respironics Deutschland

+49 8152 93 06 0

Philips Respironics France

+33 2 51 89 36 00

Philips Respironics Italy

+39 039 203 1

Philips Respironics Sweden

+46 8 120 45 900

Philips Respironics Switzerland

+41 6 27 45 17 50

Philips Respironics United Kingdom

+44 800 1300 845

www.philips.com/respironics

autoSV,AVAPS, Bi-Flex, BiPAP, OmniLab, PC Direct and Respironics are trademarks of Koninklijke Philips N.V. and its affiliates. The C-Flex mark is used under license. All rights reserved.

#### Please visit www.philips.com/omnilab



© 2013 Koninklijke Philips N.V. All rights are reserved.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

CAUTION: US federal law restricts these devices to sale by or on the order of a physician. Geyer SB 10/08/13 MCI 4105826 PN 1112817