

Appendices: Assessment of One-way Valve Efficiency in the OptiChamber Diamond VHC During Exhalation

These appendices have been produced to provide further results from the Philips Respironics white paper: “Assessment of One-way Valve Efficiency in the OptiChamber Diamond VHC During Exhalation” by Eric Lieberman. The appendices show the remaining valve leakage results for the AeroChamber Plus Z-Stat, AeroChamber FLOWSignal and ACE VHCs, as well as the full set of results for the AeroChamber Plus Flow-Vu VHC. These appendices should be read in conjunction with the main report.

Appendix A

AeroChamber Plus Z-Stat VHC results with the other breathing patterns tested:




Figure A1. Flow rate through the AeroChamber Plus Z-Stat VHC, during inhalation ■ and exhalation ■, at a rate of 15 BPM and with a tidal volume of 0.1 L. The Y axis represents airflow and the X axis represents time.




Figure A2. Flow rate through the AeroChamber Plus Z-Stat VHC, during inhalation ■ and exhalation ■, at a rate of 25 BPM and with a tidal volume of 0.5 L. The Y axis represents airflow and the X axis represents time.




Figure A3. Flow rate through the AeroChamber Plus Z-Stat VHC, during inhalation ■ and exhalation ■, at a rate of 25 BPM and with a tidal volume of 0.75 L. The Y axis represents airflow and the X axis represents time.

Appendix B

AeroChamber FLOWSignal VHC results with the other breathing patterns tested:




Figure B1. Flow rate through the AeroChamber with FLOWSignal VHC, during inhalation ■ and exhalation ■, at a rate of 15 BPM and with a tidal volume of 0.1 L. The Y axis represents airflow and the X axis represents time.




Figure B2. Flow rate through the AeroChamber with FLOWSignal VHC, during inhalation ■ and exhalation ■, at a rate of 25 BPM and with a tidal volume of 0.5 L. The Y axis represents airflow and the X axis represents time.




Figure B3. Flow rate through the AeroChamber with FLOWSignal VHC, during inhalation ■ and exhalation ■, at a rate of 25 BPM and with a tidal volume of 0.75 L. The Y axis represents airflow and the X axis represents time.

Appendix C

ACE VHC results with the other breathing patterns tested:




Figure C1. Flow rate through the ACE VHC, during inhalation ■ and exhalation ■, at a rate of 15 BPM and with a tidal volume of 0.1 L. The Y axis represents airflow and the X axis represents time.




Figure C2. Flow rate through the ACE VHC, during inhalation ■ and exhalation ■, at a rate of 25 BPM and with a tidal volume of 0.5 L. The Y axis represents airflow and the X axis represents time.




Figure C3. Flow rate through the ACE VHC, during inhalation ■ and exhalation ■, at a rate of 25 BPM and with a tidal volume of 0.75 L. The Y axis represents airflow and the X axis represents time.

Appendix D

AeroChamber Plus Flow-Vu VHC results with all the breathing patterns tested:




Figure D1. Flow rate through the AeroChamber Plus Flow-Vu VHC, during inhalation ■ and exhalation ■, at a rate of 15 BPM and with a tidal volume of 0.1 L. The Y axis represents airflow and the X axis represents time.




Figure D2. Flow rate through the AeroChamber Plus Flow-Vu VHC, during inhalation ■ and exhalation ■, at a rate of 15 BPM and with a tidal volume of 1.0 L. The Y axis represents airflow and the X axis represents time.




Figure D3. Flow rate through the AeroChamber Plus Flow-Vu VHC, during inhalation ■ and exhalation ■, at a rate of 25 BPM and with a tidal volume of 0.5 L. The Y axis represents airflow and the X axis represents time.




Figure D4. Flow rate through the AeroChamber Plus Flow-Vu VHC, during inhalation ■ and exhalation ■, at a rate of 25 BPM and with a tidal volume of 0.75 L. The Y axis represents airflow and the X axis represents time.

Appendix E

Equipment specifications for all VHCs tested:

OptiChamber Diamond

- *Manufacturer:* Respironics Respiratory Drug Delivery (UK) Ltd, a business of Philips Electronics UK Limited, Chichester, UK
- *Lot #:* 11011011
- *Valve Type:* Duckbill

AeroChamber Z-Stat

- *Manufacturer:* Monaghan Medical Corp, Plattsburgh, NY, USA
- *Lot #:* 4964-02
- *Valve Type:* Doughnut-shaped

AeroChamber Plus Flow-Vu

- *Manufacturer:* Monaghan Medical Corp, Plattsburgh, NY, USA
- *Lot #:* 203021
- *Valve Type:* Doughnut-shaped

AeroChamber with FLOWSignal

- *Manufacturer:* Monaghan Medical Corp, Plattsburgh, NY, USA
- *Lot #:* 2320-01
- *Valve Type:* Flap

ACE

- *Manufacturer:* Smiths Medical International Limited, Hythe, Kent, UK
- *Lot #:* 2422728
- *Valve Type:* Duckbill

Appendix F

Equipment specifications for the laboratory equipment used:

TSI Flow Meter

- *Manufacturer:* TSI Inc., Shoreview, MN, USA
- *Tool Calibration #:* TCN-0228
- *Calibration Due Date:* November 30, 2013
- *Serial #:* 40430749005
- *Model:* 4043 E

Pneumotach

- *Manufacturer:* Hans Rudolph Inc., Kansas City, MO, USA
- *Tool Calibration #:* TCN-0196
- *Calibration Due Date:* October 1, 2013
- *Serial #:* 381-9990
- *Model:* 4813

Flow/Volume Simulator

- *Manufacturer:* Hans Rudolph Inc., Kansas City, MO, USA
- *Serial #:* 112-024
- *Series:* 1120
- *Reference:* 113266
- *Software Version:* 4.9.4

Oscilloscope

- *Manufacturer:* Teledyne Lecroy, NY, USA
- *Tool Calibration #:* TCN-0152
- *Calibration Due Date:* October 27, 2013
- *Serial #:* LCRY0607M12644
- *Model:* WR44XI

Online interactive page can be accessed via this [\[OptiChamber Diamond Valve link\]](#), or by scanning the QR code below.



Respironics and OptiChamber are trademarks of Koninklijke Philips Electronics N.V. and its affiliates. All rights reserved.

© 2013 Koninklijke Philips Electronics N.V. All rights are reserved. Reproduction in whole or in part is prohibited without the prior written consent of the copyright holder. Part Number: 1112648. Issue 1.

