

Preferred Room Layout

12'-0" [3658mm]

Gantry €

4'-9" [1448mm]

11'-0" [3353mm]

4'-4" [1321mm]

Table €

CS

WS

TR

UPS

WT


JW

AC

CT

Tabletop to be provided by customer/contractor

Minimum Ceiling Height: 7'-0" (2134mm)



Equipment Legend				
<div> <div> <div></div> <div></div> </div> <div> <div></div> <div></div> </div> </div> <div> <div></div> <div></div> </div> <div> <div></div> <div></div> </div> <div> <div></div> <div></div> </div> <div> <div></div> <div></div> </div> <div> <div></div> <div></div> </div>				
<div> <div></div> <div></div> </div> <div> <div></div> <div></div> </div> <div> <div></div> <div></div> </div> <div> <div></div> <div></div> </div> <div> <div></div> <div></div> </div> <div> <div></div> <div></div> </div> <div> <div></div> <div></div> </div>				
Equipment Designation				
		Description	Weight lbs [kg]	Heat Load Btu/hr [W]
A	WS	CardioMD System (supine & prone)	1797[816]	2360[692]
A	WT	Patient Table		
A	AC	Acquisition Terminal (on cart)	66[30]	
A	TR	Cable Tray	- [-]	- [-]
A	CH	Calibration Source Holder Tripod	- [-]	- [-]
A	CS	Collimator Cabinet	75[35]	- [-]
G	UPS	2.2 kVA Uninterruptible Power Supply	112[51]	275[81]
G	JW	Jetstream Workstation	59[27]	1536[451]

Environmental Requirements for General Equipment Locations

The camera room HVAC system must maintain the temperature between 59°-86° F (15°-30° C) with less than 15° F (8° C) variation per hour. Humidity must be between 45% to 80%. These requirements exist on a 24 hours per day, 7 days per week basis.

Power Requirements

Input Voltage: 120 V - dedicated

Peak/Quiescent Current: 8.8 A / 5.6 A

Frequency: 50 or 60 Hz

Receptacle: NEMA L5-20 (120 V) w/ 2.2 kVA UPS

Remote Service Diagnostics

Medical Imaging equipment to be installed by Philips is equipped with a service diagnostic feature which allows for remote and on-site service diagnostics. To establish this feature, a RJ45 type Ethernet 10/100/1000 Mbit network connector must be installed. Access to customer's network via their remote access server is needed for Remote Service Network (RSN) connectivity. All costs with this feature are the responsibility of the customer.



© Koninklijke Philips Electronics N.V. 2011.
All rights reserved.

Reproduction in whole or in part is
prohibited without prior written consent of
the copyright holder.

Rev. 11.00