		E. La		
	Carlos and			
Teleman Bernan B	when the second	ALINA ALIANA ALIANA		
		241.200 E 10		
			100 200	
Ma Jaj Uzminika Kasa Kasa Majara	F 7		1000	

Validated and efficient scoring solution

Somnolyzer 24x7



sense and simplicity

Enhance your sleep study scoring

Somnolyzer 24x7 is a clinically-validated, computer-assisted sleep scoring system designed to help simplify and improve sleep study scoring productivity, accuracy, and consistency. By reducing the time sleep technicians spend on studies – typically to less than one-half of normal time – they can be free to deal with other important lab activities.

Many labs today face an increased demand for sleep diagnostic testing and a possible scarcity of sleep specialists. Implementing an automated scoring solution, like Somnolyzer 24x7, can help deal with these challenges. Automated scoring can help you enhance your lab efficiency, forecasting, and demand. Designed by sleep experts who understand your needs, Somnolyzer 24x7 offers many benefits for lab managers and technicians. It also can help enhance the relationship between labs and referring physicians and their patients who want quick study results.

UlireSh	na jana po jak (m. et et	anna an	S ng		AL PULL	
	and the first matter of the first of the fir	and the second and the	Million and	\$		
	man which when the	Il mlitter	auwww.wool	-	9	
	man Muchard	when the when the state of the	undertaille marine	2	No.	
	All draw and from		he h			
		1 (1 , 11), 1), 1 7, 1 7, 1 , 1		distant.	-	
		Deal	Same			

Provides a smart solution to the challenges of scoring sleep studies

Somnolyzer 24x7 has been proven to most closely approximate human scoring and was tested on the largest normative database, as well as on thousands of patient studies.¹

Some technicians may find they only need as little as 10-20 minutes for structured quality expert review. This increased efficiency may enable them to focus their expertise on the areas requiring refinement instead of routine, monotonous scoring.

Somnolyzer 24x7 is a reliable and validated, computerassisted sleep scoring system that:

• Includes an expert review process

Overall performance (staging)


- Applies the same scoring rules as any human expert
- Follows R&K and current AASM standards
- · Detects all relevant visual features (i.e., spindles, eye movements, delta waves) just as in manual scoring


- · Enables you to view reports in your native system without having to change them (dependent upon the sleep system you use)
- Ensures that your scored results are of consistent quality by diminishing interscorer variability

Somnolyzer 24x7 provides automated analysis of:

- Sleep staging
- Arousal analysis
- Respiratory analysis
- Leg movement analysis

Somnolyzer 24x7 is designed to help your scoring studies be more efficient, consistent, and maintain a high level of inter-scorer agreement. With Somnolyzer 24x7 in place, your staff should be able to provide referral sources and patients with quicker results and handle any possible demands for increased patient throughput.





Note: the above graphs are adapted from the studies referenced in the footnotes.

Widely recognized for its utility

Somnolyzer 24x7 adheres to the "AASM Manual for the Scoring of Sleep and Associated Events" 2007 guidelines.

- It was listed as the top Digital Sleep Analysis solution by the AASM in their evidence based review¹ of Sleep Stage validation (as detailed in a multi-center study with the largest number of subjects)
- It has been used and validated in more than 20 drug and device clinical trials (insomnia, narcolepsy, etc.)^{4,5,6}


Enhance your lab operation

your services

Somnolyzer 24x7 provides lab managers with the ability to enhance:

Efficiency – enables a more efficient use of staff resources Productivity – provides a significant increase in throughput with the same resources and helps reduce backlogs Planning – easier capacity and financial planning with more

predictable staffing requirements Growth – by increasing scoring efficiency you can be better positioned to address increased demand or expand



The Polysomnographic Feature File (PFF) helps clinicians easily identify signal deficiencies so they can focus their expertise where it is needed most. Trends within the file depict the distribution of key features in the EEG, EOG or EMG. The events (respiratory, leg movements, arousals) review is based upon an all-night overview of events and on a quick review based on a five-minute window of both signals and events. Looking at single epochs or single events is only necessary in clearly identified cases, based on this global picture.

Satisfy referral sources

Somnolyzer 24x7 enables you to provide many benefits to referring physicians:

Reliability – achieve close to 100% inter-scorer agreement following expert review^{2,3}

Continuity – no need to adapt to a new layout; depending upon your diagnostic system, reports can look the same as before

Turnaround – reports are received much more quickly, sometimes even the same day

Standardization – patient results are compliant with the latest professional guidelines

Highly efficient, accurate, and consistent

Just as with human scoring, Somnolyzer 24x7 is most accurate when the signals are of high quality and compliant with the standard.

With its computer-supported solution followed by an expert review, Somnolyzer 24x7 helps enable the greatest degree of accuracy and consistency without any systematic deviations that could have an influence on subsequent diagnosis.

Works with other PSG systems

Even if your diagnostic system is not fully supported for reintegration of Somnolyzer 24x7 results, our comprehensive sleep report – the "Spot Report" – contains all major variables (from AHI to Total Sleep Time).

All variables are compared to the world's largest normative database on healthy sleep. Depending on the age and gender of the subjects a sleep profile of the most important parameters are compared to the norm. Significant deviations are clearly indicated and marked, and the report is generated in PDF format.


Integrated work flow

Somnolyzer 24x7 can directly interface with Alice Sleepware so that the user does not have to change work flow. Somnolyzer 24x7 results also can be reintegrated into a number of other leading diagnostic systems.

Home testing capability

Somnolyzer 24x7 can provide scoring on ambulatory or home testing studies for:

- Full portable PSG (Level 2; ambulatory PSG)
- Cardio-respiratory testing (Level 3)
- One- or two-channel devices (Level 4)



Part number	Description
1076888	SOMLZ PSG Analysis (volume-based pricing)
1076892	SOMLZ Onsite Server Install/Training
1076893	SOMLZ PR Hosted Server Install/Remote Training
1076890	SOMLZ Expert Review Service PSG
1076891	SOMLZ Non-Contract Full Service

For home sleep testing (ask your sales representative regarding availability)		
Part number	Description	
1076889	SOMLZ HST Analysis (volume-based pricing)	
1078155	SOMLZ Expert Review Service HST	

Philips Healthcare is part of Royal Philips Electronics

How to reach us www.philips.com/healthcare healthcare@philips.com

Asia +49 7031 463 2254

Europe, Middle East, Africa +49 7031 463 2254

Latin America +55 11 2125 0744

North America +1 425 487 7000 800 285 5585 (toll free, US only) Philips Respironics 1010 Murry Ridge Lane Murrysville, PA 15668

Customer Service +1 724 387 4000 800 345 6443 (toll free, US only)

Philips Respironics International Headquarters +33 1 47 28 30 82

Philips Respironics Asia Pacific +65 6882 5282

Philips Respironics Australia +61 (2) 9666 4444

Philips Respironics China +86 021 24127311 Philips Respironics Deutschland +49 8152 93 06 0

Philips Respironics France +33 2 51 89 36 00

Philips Respironics Italy +39 039 203 1

Philips Respironics Sweden +46 8 120 45 900

Philips Respironics Switzerland +41 6 27 45 17 50

Philips Respironics United Kingdom +44 800 1300 845

www.philips.com/respironics

- ¹ Penzel T, Hirshkowitz M, Harsh J, et al. Digital analysis and technical specifications, Journal of Clinical Sleep Medicine 2007;3(2)109-120.
- ² An E-Health Solution for Automatic Sleep Classification According to Rechtschaffen and Kales: Validation Study of the Somnolyzer 24x7 Utilizing the Siesta Database; Neuropsychobiology 2005;51:115–133.
- ³ Computer-assisted Sleep Classification According to the Standard of the American Academy of Sleep Medicine: Validation Study of the AASM Version of the Somnolyzer 24x7; Neuropsychobiology 2010;62:250–264
- ⁴ Svetnik V, Ma J, Soper K et al. Evaluation of automatic and semi-automatic scoring of polysomnographic recordings from a clinical trial using Zolpidem in the treatment of insomia, Sleep, 30(11):1562-74.
- ⁵ Dingemanse J, Dorffner G, Hajak G, et al. Proof-of-concept study in primary insomnia patients with almorexant (ACT-078573), a dual orexin receptor antagonist, Sleep and Biological Rhythms, Worldsleep07 supplement, 2007.
- ⁶ Mayleben D, Safirstein B, Rajachandran L, et al. Aplindore, a partial dopamine agonist, reduces periodic limb movements in patients with restless legs syndome. Sleep 2009; 32, abstract supplement, A304.

Respironics, Alice, Sleepware, and Somnolyzer are trademarks of Respironics Inc., and its affiliates. All rights reserved.

Please visit www.philips.com

PHILIPS

©2010 Koninklijke Philips Electronics N.V. All rights are reserved.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

CAUTION: US federal law restricts these devices to sale by or on the order of a physician.

Hstetter SB 09/28/10 MCI 4103543 PN 1078246