
# Determining the Influence of Washing on the Aerosol Performance of an Anti-Static Valved Holding Chamber

D. von Hollen<sup>1</sup>, L. Slator<sup>2</sup>, K. Nikander<sup>1</sup>, R. H. M. Hatley<sup>2</sup>. <sup>1</sup> Philips Respironics, Respironics New Jersey, Inc., Parsippany, New Jersey, USA. <sup>2</sup> Respironics Respiratory Drug Delivery (UK) Ltd, Chichester, West Sussex, UK.



effective aerosol delivery to respiratory patients. This study compared the aerosol characteristics of an unwashed ('outof-the-box') and prewashed preproduction Diamond VHC using three pMDI drug formulations.


Figure 4. Mean fine particle dose (amount of drug in NGI  $\leq$  4.7 µm) from the unwashed Diamond VHCs and prewashed

Fluticasone propionate


Ipratropium bromide


Anti-static VHC

#### **Conventional VHC**



If there is a build up of static charge on the VHC walls, upon actuation the aerosol particles with opposite charges will be attracted to the VHC walls, leading to a rapid reduction in the dose available for inhalation.





VHC keeps aerosol in chamber longer giving time for patient to inhale complete dose of drug.

Increased convience for user - washing requirement for hygiene.

User should wash with ionic detergent frequently washing reduces static charge.

Without regular and frequent washing the dose

available for inhalation reduces rapidly over time.

Diamond VHCs. Error bars denote standard deviation about the mean.

Albuterol sulfate




Figure 5. Mean MMAD (mass median aerodynamic diameter) from the unwashed Diamond VHCs and prewashed Diamond VHCs. Error bars denote standard deviation about the mean.


## Discussion

0

Both the fine particle doses and MMADs for the three pMDI drug formulations were similar from unwashed and prewashed Diamond VHCs. The practical benefit of these results is that they show that users and care givers can use a new Diamond VHC for treatment as soon as it is removed from the packing without having to wash the device first, as is the case with conventional VHCs.

#### Conclusions

## Method



#### Figure 3. Experimental test method.

After each test, the induction port, back-up filter, NGI cups and VHCs were processed using the relevant HPLC assay diluent. CITDAS V3.10 software was used to generate the aerosol characteristics data. The fine particle dose (amount of drug in NGI  $\leq$  4.7 µm), and Mass Median Aerodynamic Diameter (MMAD) were calculated.

- There were minimal differences in the aerosol characteristics between the unwashed and prewashed pMDI VHC combinations with all of the drugs tested.
- The results indicate that the aerosol performance of the Diamond VHC is the same when used in either unwashed or prewashed states.

## References

- 1) Virchow J.C., Crompton G.K., Dal Negro R., Pederson S., Magnan A., Seidenberg J., Barnes P.J. Importance of inhaler devices in the management of airway disease. Respir Med. 2008;102:10-19.
- 2) Piérart F., Wildhaber J.H., Vrancken I., Devadason S.G. Le Souëf. Washing plastic spacers in household detergent reduces electrostatic charge and greatly improves delivery. Eur Respir J. 1999;13:673-678
- 3) Kwok P.C.L., Chan H. Electrostatic charges of metered dose inhaler aerosols sampled from plastic spacers. In: Dalby R.N., Byron P.R., Peart J., Suman J.D., Farr S.J. (eds) Respiratory Drug Delivery X. 2006. Virginia Commonwealth University, Richmond, VA:901-903.
- 4) Kwok P.C.L., Glover W., Chan H. Electrostatic charge characteristics of aerosols produced from metered dose inhalers. J Pharm Sci. 2005; 95(12):2789-2799.

Acknowledgements: Drafting and editorial assistance was provided by N. Smith and L. Pearce of PS5 Consultants Ltd.

Presented at the Association of Asthma Educators Annual Conference, 22<sup>nd</sup>-24<sup>th</sup> July, 2011, Denver, Colorado, USA.





