





Dameca A/S Philips Anesthesia Care	Product Recycling Passport	Page 1 of 4
---------------------------------------	-----------------------------------	-------------

Product name:	IntelliSave AX700	
Identification code(s)	10623-00	
Total weight (in Kg)	Approx. 150 kg (depending on configuration)	
Producer/ Manufacturer	Name company:	Dameca A/S
	Address:	Islevdalvej 211; Roedovre
	Zip code:	2610
	Country:	Denmark

Recycle Info	Items:	Location
Special attention 	WARNING High Voltage - Voltages dangerous to life are present in the machine when it is connected to the mains power supply. Do not perform any disassembly procedures with power applied to the machine. Failure to adhere to this warning could cause serious injury or death.	Figure 1
	Before starting any maintenance or repair, release all gas pressures; first N ₂ O, then O ₂ and Air	Figure 2
	WARNING When removing and mounting the battery terminals and while they are unmounted, take care not to short-circuit the terminals with the metal cover	
 Fluids / Gases	Items: Absorber material	Figure 3
Batteries  To be Removed	Type: 2 12 V Lead Acid batteries (EGM batteries)	Figure 4
	20-Cell NiMH, 24V, 1.5Ah, (ventilator battery)	Figure 4
	Polycarbonmonofluride lithium battery (PCGM battery)	Figure 4
	Manganese dioxide lithium battery (Log PC battery)	Figure 4
Hazardous  To be Removed	Substances:	Location
	PCB boards contain leaded solder (<i>material data is needed</i>):	Figure 4
	Lead-acid batteries weight in total is 5 kg	Figure 4

Title: Product Recycling Passport for AX700
 Recycling Passport Number: 10623 RCP Rev: 50

Copies are uncontrolled

Show locations of materials mentioned on the previous sheet.

Figure 1:
Location and How to disconnect power



Figure 2: Location and how to release pressure



Push to release pressure from walloutlet (AGA).



Turn nut counterclockwise to release from machine inlet (NIST)

Figure 3:




Figure 4:

