
Reusable SpO₂ T-Series Sensors

Instructions for Use

M1191T

Adult Finger Sensor

M1192T

Pediatric Finger Sensor

M1193T

Neonate Hand, Foot Sensor

PHILIPS

About This Edition

Notice

Copyright © 2007

Koninklijke Philips Electronics N.V.
All rights are reserved.

Printing History

Current Edition/Revision:

Printed in USA January 2007

Previous Editions/Revisions:

Printed in Germany 07/2005

Printed in Germany 01/2004

Contents

Instructions for Use.....	1	Návod na použitie	127
Manuel d'utilisation	10	Navodila za uporabo	136
Gebrauchsanweisung	19	Instrucțiuni de utilizare	145
Gebruiksaanwijzing	28	Инструкция за ползване	154
Instrucciones de uso	37	Kasutusjuhend	163
Instruções de Uso	46	Naudojimo instrukcija.....	172
Istruzioni per l'uso	55	Lietošanas pamācība	181
Bruksanvisning	64	Инструкция по эксплуатации. 190	
Bruksanvisning	73	Οδηγίες Χρήσης.....	199
Käyttöopas	82	Kullanım Talimatları.....	208
Brugerhåndbog.....	91	ユーザーズ・ガイド	217
Használati útmutató	100	使用説明.....	226
Instrukcja obsługi.....	109	使用说明书.....	235
Návod na použití	118		

USA Restriction

United States law restricts this device to sale by or on the order of a physician.

Product Overview

Intended Use

Philips Reusable SpO₂ T-Series sensors provide continuous, non-invasive measurement of arterial oxygen saturation (provides pulse rate monitoring and plethysmograph wave) to compatible SpO₂ patient monitoring devices.

The sensors are suitable for use only by appropriately trained clinicians such as nurses, physicians, and EMS personnel for the treatment of patients in and out of a hospital, as well as during transport. All clinicians shall follow the current standard use model when applying a SpO₂ sensor to a patient; i.e. apply the M1191T Adult or M1192T Pediatric sensor to a patient's finger; apply the M1193T Neonate sensor to the patient's hand or foot (see table below).

T-Series Sensor/Monitor Compatibility

T-Series sensors are compatible with many *validated* Philips/Agilent/HP and OEM (Original Equipment Manufacturer) monitors. Reference the *Sensor & Instrument Compatibility Sheet* included with your T-Series Sensor for a comprehensive list of compatible monitors.

Reusable T-Series Sensors	Model	Application	Cable	Connector
	M1191T Sensor	Adult Any finger except thumb	45 cm	 9-pin D-sub Connector
	M1192T Sensor	Pediatric Any finger except thumb	45 cm	
	M1193T Sensor	Neonatal Foot or Hand	90 cm	

Accessories Supplied with Sensors

- **Sensor & Instrument Compatibility Sheet:** Lists all instruments that are compatible with the reusable T-Series Sensors.
- **Instructions for Use:** Includes T-Series Sensor product descriptions, patient applications, Warnings against improper use, cleaning and disinfection procedures, and sensor specifications.

Accessories Available for Separate Purchase

- **M1900B Adapter Cable:** 9-pin (D-Sub female) to 12-pin (round male) *adapter* cable. Connects any T-Series Sensor 9-pin connector to any compatible Philips/Agilent/HP instrument that has a 12-pin SpO₂ connection. Extends sensor cable length by 3m.
- **M1943A Adapter Cable:** 9-pin (D-sub female) to 8-pin (round male) *adapter* cable. Connects any T-Series Sensor 9-pin connector to any compatible Philips/Agilent/HP instrument that has an 8-pin SpO₂ connection. Extends sensor cable length by 1.1m.
- **M1943AL Adapter Cable:** A longer (3m) version of the M1943A *adapter* cable listed above.

Definitions of Product Symbols

					
Adult Patient	Pediatric Patient	Neonate Patient	Neonate Foot	Storage Temperature Limitations	Storage Humidity Limitations
			R _X only		
Adult Finger	Pediatric Finger	Neonate Hand	Prescription Use Only		

Warnings

- These sensors are for use only with compatible monitoring instruments. Before using, reference your instrument's Instructions for Use, or the *Sensor & Instrument Compatibility Sheet* shipped with each sensor, to verify instrument compatibility. Be sure to use sensors only with validated instruments, otherwise patient injury may result.
- Unapplied sensors might cause readings. To avoid misdiagnosis, verify sensor is applied to patient correctly.
- Pulse oximetry measurements are statistically distributed. Two-thirds of all pulse oximetry measurements can be expected to fall within the stated accuracy (refer to *Specifications* later in this IFU, or the *Sensor & Instrument Compatibility Sheet* for stated sensor accuracy).
- Connect the sensor only to the SpO₂ connection or the SpO₂ adapter cable of the oximeter.
- Do not reuse the sensor on a *different* patient until it has been disinfected. The sensor can be reused on the same patient for the patient's entire stay.
- At elevated ambient temperatures, patient skin could be severely burned after prolonged sensor application at sites that are not well perfused. To prevent this condition, be sure to check patient application sites frequently. All listed sensors operate without risk of exceeding 41° C on the skin if the initial skin temperature does not exceed 35° C.
- Be sure to apply the sensor correctly at the preferred application site, following instructions provided later in this document (refer to *Applying the Sensor*). Failure to do so can cause inaccurate measurements.
- To avoid venous pulsation, obstructed circulation, pressure marks, pressure necrosis, artifacts and inaccurate measurements, make sure you are using a sensor of the correct size and that the sensor is not too tight. If the sensor is too tight, because the application site is too large or becomes too large due to edema, excessive pressure may be applied. This can result in venous congestion distal from the application site, leading to interstitial edema and tissue ischemia.

Warnings (continued)

- If the sensor is applied too loose, it may fall off or compromise proper alignment of sensor optics and result in inaccurate readings.
- Where possible, the application site for the sensor should be an extremity free of arterial catheters, blood pressure cuffs, or intravascular infusion lines.
- Avoid application sites subject to excessive motion. Try to keep the patient still, or move the sensor to a site with less motion.
- Dysfunctional hemoglobin or intravascular dyes can cause inaccurate measurements.
- Make sure the application site for the sensor is not deeply pigmented, or deeply colored (for example: nail polish, artificial nails, dye or pigmented cream may cause inaccurate measurements). In any of these cases reposition the sensor or choose an alternative sensor for use at a different site.
- Cover the sensor with opaque material under conditions of strong or excessive light (infrared lamps, OR lamps, photo therapy). Failure to do so can result in inaccurate measurements.
- Protect the connector from contact with any liquid.
- Inspect the sensor application site every 2 to 3 hours to ensure skin integrity, correct optical alignment, and circulation distal to the sensor site. Skin irritation or ulceration may occur if the sensor is attached to one location for too long. Move the sensor application site every 4 hours, or more often if circulation or skin integrity is compromised. If the light source is not directly opposite the light detector, reapply the sensor, or choose an alternative sensor for use at a different site.
- Do not use a sensor during MRI scanning. This may cause burns or inaccurate measurements.

Before Applying the Sensor

Before applying any T-Series Sensor to a patient, be sure to read and understand all Warnings listed in your SpO₂ monitoring instrument IFU, as well as the Warnings in this sensor IFU. Always apply an appropriate size sensor at the proper patient application area (finger for adult and pediatric patients; hand or foot for neonatal patients).

Inspect Sensor for Damage

- Inspect the sensor outside and inside. To check the inside, gently open the sensor cavity and check for splits on or next to the transparent silicone that covers the optical elements. Be sure there are no blisters on the silicone, and no liquid leaking from sensor optics.
- Any sensor showing signs of damage or alteration must not be used for further patient monitoring; instead, dispose of it using proper disposal procedures (see below).

Sensor Disposal

Dispose Deteriorated Sensors

Any sensor showing signs of physical or electrical deterioration or failure should be disinfected, decontaminated, and disposed of according to local laws regarding the disposal of hospital waste.

Change Application Site Periodically

Move the sensor application site every 4 hours, or more often if circulation or skin integrity is compromised.

Applying the Finger Sensor

M1191T Adult Finger Sensor

For patients weighing more than 50kg.

M1192T Pediatric Finger Sensor

For patients weighing between 15kg and 50kg.

*adult application shown
with supplied wristband*



Step	Applying the M1191T (Adult) or M1192T (Pediatric) Sensor
1	Select appropriate sensor for patient size (defined above).
2	Position sensor over patient's finger, being sure the sensor cable is positioned ON TOP of the finger/hand as shown above.
3	Patient fingertip should touch but not protrude from end of sensor. If necessary, trim the fingernail to correctly position the sensor.
4	If applying a M1191T Adult finger sensor, secure the cable to BACK OF HAND with supplied M1627A wristband (wristband is supplied only with M1191T Adult finger sensors.)
5	Plug sensor into instrument (or into adapter cable if needed).
6	Inspect and change the sensor application site periodically. You may reapply the sensor to the <i>same</i> patient by repeating this procedure. Clean and disinfect the sensor before applying it to a <i>different</i> patient.

Applying the Hand/Foot Sensor

M1193T Neonatal Hand/Foot Sensor

For patients weighing 1kg-4kg.

Step	Applying the M1193T (Neonate) Sensor
1	Position the sensor over the hand or foot, with Sensor Optics aligned opposite each other. Hold the sensor in position while stretching the Strap slightly, not more than 2.5cm (1.0"). 
2	Press the stretched Strap against the Slot, thread it into the Latch, and pull for a secure fit. If long enough, thread it through the Second Latch as well. 
3	Plug the sensor cable/connector into the monitoring instrument (or the instrument's adapter cable if used).
4	Inspect and change the sensor application site periodically. You may reapply the sensor to the <i>same</i> patient by repeating this procedure. Clean and disinfect the sensor before applying it to a <i>different</i> patient.

Sensor Cleaning and Low-Level Disinfection

The reusable sensors should be cleaned and disinfected, but never sterilized. Follow the procedure outlined below.

Warnings

- Use only validated cleaning agents and disinfectants listed below; do not use any others. Failure to do so may damage the sensor or its connecting wires and shorten product lifetime or cause safety hazards.
- Select disinfectants carefully as some have very similar names but completely different compositions.
- Do not immerse the sensor connector in any of the cleaning solutions, disinfectant, or other liquid (only the sensor and cable housing may be immersed, not the connector).
- Do not soak sensors in disinfectants longer than specified by the disinfectant manufacturer.
- Do not sterilize the sensors.

Validated Cleaning Agents

- Mild Detergent
- Salt Solution (1%)

Validated Disinfectants

Metricide [®] 28	Cidex [®] Formula 7	Kohrsolin [®] (2%)
Metricide [®] Plus 30	Cidex [®] OPA	Mucapur [®] -CD (1%)
Terralin [®] Liquid	Cidex [®] Plus	Isopropanol (70%) or Isopropanol Wipe (70%)
Incidin [®] Liquid	Omnicide [®] 28	

Step	Cleaning and Low-Level Disinfection
1	Clean sensor following instructions supplied with cleaning agent.
2	Disinfect sensor following instructions supplied with disinfectant.
3	Rinse sensor in water, wiping it dry with a clean cloth, leaving to dry completely. If you notice any deterioration or damage to sensor or cable, dispose the sensor immediately.

Specifications

Sensor Accuracy

- M1191T (Adult) and M1192T (Pediatric) sensor accuracy is $\pm 3\%$ SpO₂ in the range from 70% to 100% SpO₂.
- M1193T (neonate) sensor accuracy is $\pm 3\%$ SpO₂ in the range from 70% to 100% SpO₂. When sensors are used on neonatal subjects as recommended, the literature predicts the accuracy range should be increased by $\pm 1\%$ to account for affect on oximeter measurements due to fetal hemoglobin in neonatal blood.

For more information on SpO₂ and Pulse Rate accuracy specifications, reference your monitoring instrument's Instructions for Use, or the *Sensor & Instrument Compatibility Sheet* included with your sensor.

Light Emitting Diode Wavelength Ranges

Wavelength ranges for the light emitting diodes used in these sensors are within 600nm-1000nm, with an optical output power of less than 15mW. Knowing the wavelength range can be useful to clinicians performing photodynamic therapy.

Measurement Validation

The SpO₂ accuracy has been validated in human studies against arterial blood sample reference measured with a CO-oximeter. In a controlled desaturation study, healthy adult volunteers with saturation levels between 70% and 100% SaO₂ were studied. The population characteristics for those studies were:

- Approximately 50% female and 50% male ranging in age from 18-45
- Skin tone: from light to black

Because pulse oximeter equipment measurements are statistically distributed, only approximately 2/3 of pulse oximeter equipment measurements can be expected to fall within the $\pm$ Arms value measured by a CO-oximeter. Functional testers, such as an SpO₂ simulator cannot be used to assess the accuracy of pulse oximeter sensors.

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For more information, please call your local Philips sales office listed in your telephone directory or a Philips regional office listed below for the location of your nearest sales office.

United States:

Philips Medical Systems
North America Corporation
3000 Minuteman Road
Andover, MA 01810-1099
(800) 225-0230

Canada:

Philips Medical Systems Canada
281 Hillmount Road
Markham, ON
L6C 2S3
(800) 291-6743

Latin America Headquarters:

Philips Medical Systems
1550 Sawgrass Corporate Parkway #300
Sunrise, FL 33323 USA
Tel: (954) 835-2600
Fax: (954) 835-2626

Manufacturing Address:

Philips Medizin Systeme
Böblingen GmbH
Hewlett-Packard Str. 2
71034 Böblingen, Germany

Europe, Middle East and Africa:

Philips Medizin Systeme
Böblingen GmbH
Hewlett-Packard Str. 2
71034 Böblingen, Germany
Fax: (+49) 7031 463 1552

Asia Pacific Headquarters:

Philips Medical Systems
30/F Hopewell Centre
17 Kennedy Road
Wanchai
Hong Kong
Tel: (852) 2821 5888
Fax: (852) 2527 6727

<http://shop.medical.philips.com>

 **0366** This product complies with the requirements of the Council Directive 93/42/EEC June 1993 (Medical Device Directive).

 Contains no Latex

Subject to modification

© 2007

Koninklijke Philips Electronics N.V.
Printed in USA January 2007
Philips Part Number 453564052301