


**PHILIPS
MEDICAL
SYSTEMS**

WEEE Directive: Recycling Passport


Document # M4993-90501
Page 1 of 4
2005 July: Rev A

Product name:	PageWriter Trim Cardiographs (Trim II, III, Rx)	
Identification code(s)	PW Trim III: 860286 PW Trim II: 860288 PW Trim Rx: 860297	
Total weight (in Kg)	7.4 Kg (Cardiograph w/battery and PIM), 1.2 Kg (Battery)	
Producer/ Manufacturer	Name company:	Philips Medical Systems
	Address:	3000 Minuteman Road, Andover Ma
	Zip code:	01810
	Country:	USA

Recycle Info	Items:	Location
Special attention 	First Steps: Disconnect power cord and remove battery before starting disassembly of the unit. Special Tools Needed: Torx 10 Driver to remove LCD display, Torx 6 for Keyboard and Torx 15 for Lower case. PW Trim Repair manual 453564010431 provides additional information for disassembly.	N/A
Fluids / Gases 	Items: None	N/A
Batteries 	Type: Lead Acid Battery Philips Part # 989803130151	See page 4
To be Removed		
Hazardous 	Substances: HG in Liquid Crystal Display (LCD) with Cold cathode fluorescent lamp (CCFL) light source. Cardiograph has 6 circuit boards Patient Interface Module (PIM) has 2 circuit boards	See pages 2 Several locations inside cardiograph and PIM
To be Removed		

Show locations of materials mentioned on the previous sheet. (Picture information)

Figure 1-3 PageWriter Trim II, III and Rx Cardiograph and Cart (Front View)



- A** Patient Interface Module (PIM)
- B** Printer paper/report storage slot
- C** Storage Shelf
- D** AC Power Cord
- E** PIM Leads


- F** Wheel Brake
- G** Control Panel
- H** Keyboard
- I** Trim Knob
- J** Printer Paper Drawer



Removing the bezel




After removal of the bezel



Use Torx 10 driver to remove screws at top of display.

Figure 1-4 PageWriter Trim II, III, Rx, and CDIC Cardiograph (Rear View)



- | | |
|-------------------------------|---|
| K Reset button | Q Battery door |
| L PC card eject button | R PIM connector |
| M PC card slot | S SmartCard Reader or USB memory stick connector |
| N AC power connector | T Barcode reader or magnetic card reader connector |
| O Equipotential post | U Serial connector (not supported) |
| P Fuse door | |

