Reusable SpO₂ Sensors

Instructions for Use

M1191B and M1191BL

Adult Finger Glove

M1192A

Pediatric Finger Glove

M1193A

Neonate Hand/Foot Wrap

M1194A

Adult/Pediatric Ear Clip

M1195A

Infant Finger/Toe Glove

PHILIPS

About This Edition

Notice

Copyright © 2009

Koninklijke Philips Electronics N.V.
All rights are reserved.

Printing History

Current Updated Edition/Revision:

March 2009
Part Number 453564154671

Previous Editions/Revisions:

April 2008
Part Number 453564095901

December 2006
Part Number 453564052251

Contents

Instructions for Use.....	1	Návod na použitie	127
Manuel d'utilisation	10	Navodila za uporabo	136
Gebrauchsanweisung	19	Instrucțiuni de folosire	145
Gebruiksaanwijzing	28	Указания за употреба.....	154
Instrucciones de uso	37	Kasutusjuhend.....	163
Instruções de Uso	46	Naudojimo instrukcija.....	172
Istruzioni d'uso	55	Lietošanas instrukcijas	181
Brøkerhåndbok.....	64	Инструкция по эксплуатации ...	190
Bruksanvisning	73	Οδηγίες Χρήσης.....	199
Käyttöopas	82	Kullanım Talimatları.....	208
Brugsanvisning	91	ユーザーズ・ガイド.....	217
Használati útmutató	100	使用説明.....	226
Instrukcja obsługi.....	109	使用说明书.....	235
Návod k obsluze.....	118		

United States law restricts this device to sale by or on the order of a physician.

Product Overview

Philips A-Series and B-Series reusable SpO₂ sensors are available in three different styles and can be used to monitor a broad range of patient types as indicated in the *Sensor Applications* table below. Refer to *Specifications* later in this IFU for sensor accuracy.

Intended Use

These SpO₂ sensors are for multi-patient use, when continuous non-invasive arterial oxygen saturation and pulse rate monitoring are required.

Instrument Compatibility

All A-Series sensors (model numbers ending with A) are for use with any Philips/Agilent/HP instrument that lists those sensors as an accessory. The B-Series sensors (model numbers ending with B or BL) are for use with any Philips/Agilent/HP instrument that lists them, or the obsoleted M1191A and M1191AL sensors, as compatible accessories; i.e. the M1191B sensor replaces the obsoleted M1191A; the M1191BL sensor replaces the obsoleted M1191AL.

Sensor Applications

Sensor Model Number	Sensor Style	Adult		Pediatric		Infant		Neonate	
		Finger	Ear	Finger	Ear	² Finger	Toe	² Hand	Foot
M1191B	Glove Style	✓							
¹ M1191BL	Glove Style	✓							
M1192A	Glove Style			✓					
M1193A	Wrap Style							✓	✓
M1194A	Ear-Clip Style		✓		✓				
M1195A	Glove Style					✓	✓		

¹M1191BL sensors have a Longer 3m (versus standard 2m) attached cable.

²Preferred application sites for M1193A and M1195A sensors.

Accessories Supplied with Sensors

- **M1627A Adult Wristband:** *Supplied only with adult glove style sensors.* Used to secure the sensor cable to top of adult patient's hand and wrist.

Accessories Available for Separate Purchase

- **M1940A Adapter Cable:** 8-pin (round female) to 12-pin (round male) *adapter* cable. Adapts 8-pin sensor connector to 12-pin instrument jack. Extends sensor cable length by 2m. DO NOT USE with M1191BL sensors.
- **M1941A Extension Cable:** 8-pin (round female) to 8-pin (round male) *extension* cable. Extends sensor cable length by 2m. DO NOT USE with M1191BL sensors.

Definition of Product Symbols

 Adult	 Pediatric	 Adult Finger	 Infant/Neo Hand	 Adult/Pediatric Ear
 Infant	 Neonate (Neo)	 Pediatric/Infant Finger	 Infant Toe	 Infant/Neo Foot
 Caution Consult Accompanying Documents	R_x only Prescription Use Only	 Storage Temp	 Contains No Latex	

Warnings

- These sensors are for use only with Philips/Agilent/HP instruments. Before using, verify that the combination of instrument/sensor is specified in the instrument's user documentation (e.g. Instructions for Use); otherwise patient injury may result.
- Pulse oximetry measurements are statistically distributed. Two-thirds of all pulse oximetry measurements can be expected to fall within the stated accuracy (refer to *Specifications* later in this IFU for stated sensor accuracy).
- Connect the sensor only to the SpO₂ connection or the SpO₂ adapter cable of the oximeter.
- Do not reuse the sensor on a different patient until it has been disinfected. The sensor can be reused on the same patient for the patient's entire stay.
- At elevated ambient temperatures, patient skin could be severely burned after prolonged sensor application at sites that are not well perfused. To prevent this condition, be sure to check patient application sites frequently. All listed sensors operate without risk of exceeding 41° C on the skin if the initial skin temperature does not exceed 35° C.
- Be sure to apply the sensor to the patient at the preferred or alternative application site, while following application instructions provided later in this document. Failure to do so can cause inaccurate measurements.
- To avoid venous pulsation, obstructed circulation, pressure marks, pressure necrosis, artifacts and inaccurate measurements, make sure you are using a sensor of the correct size and that the sensor is not too tight. If the sensor is too tight, because the application site is too large or becomes too large due to edema, the excessive pressure applied may result in venous congestion distal from the application site, leading to interstitial edema and tissue ischemia.
- If the sensor is applied too loosely it may fall off or compromise proper alignment of sensor optics and result in inaccurate readings.

Warnings (cont.)

- Where possible, the application site for the sensor should be an extremity free of arterial catheters, blood pressure cuffs, or intravascular infusion lines.
- Avoid sites subject to excessive motion. Try to keep the patient still, or move the sensor to a site with less motion.
- Dysfunctional hemoglobin or intravascular dyes can cause inaccurate measurements.
- Make sure the application site for the sensor is not deeply pigmented, or deeply colored. For example, nail polish, artificial nails, dye or pigmented cream may cause inaccurate measurements. In any of these cases reposition the sensor or choose an alternative sensor for use at a different site.
- Cover the sensor with opaque material under conditions of strong or excessive light (infrared lamps, OR lamps, photo therapy). Failure to do so can result in inaccurate measurements.
- Protect the connector from contact with any liquid.
- Inspect the sensor application site every 2 to 3 hours to ensure skin integrity, correct optical alignment, and circulation distal to the sensor site. Skin irritation or ulceration may occur if the sensor is attached to one location for too long. Change the sensor application site every 4 hours, or more often if circulation or skin integrity is compromised. If the light source is not directly opposite the light detector, reapply the sensor, or choose an alternative sensor for use at a different site.
- Do not use a sensor during MRI scanning. This may cause burns or inaccurate measurements.
- **M1191BL only:** Do not use the M1191BL sensor with an extension or adapter cable.
- **M1195A only:** Infant patients are prone to movement; therefore, when applying the M1195A sensor to an infant, secure the sensor cable with tape to prevent the sensor from coming off the infant's finger. If necessary, position the sensor on a different finger or an alternate site (toe).

Before Applying the Sensor

Be sure to read and understand all Warnings listed in your SpO₂ monitoring instrument IFU, as well as all Warnings described in this sensor IFU. Use sensors only with validated instruments and at a recommended patient application area.

Inspect Sensor for Damage

- Before reusing, inspect the sensor outside and inside. To check the inside of a glove style sensor, gently open the finger cavity and check for splits on or next to the transparent silicone that covers the optical elements. Be sure there are no blisters on the silicone and no liquid leaking from sensor optics.
- Any sensor showing signs of damage or alteration must not be used for further patient monitoring; instead, dispose of it using proper disposal procedures (see below).

Sensor Disposal

Disposal of Deteriorated Sensors

Any sensor showing signs of physical or electrical deterioration or failure should be disinfected, decontaminated, and disposed of according to local laws regarding the disposal of hospital waste.

Change Application Site Periodically

Move the sensor application site every 4 hours, or more often if circulation or skin integrity is compromised.

Applying Glove Style Sensors

Adult Patients (>50kg)

Apply M1191B or M1191BL Sensor to any finger except thumb.

*adult application
shown*

Pediatric Patients (15-50kg)

Apply M1192A Sensor to any finger except thumb.

Infant Patients (4-15kg)

Apply M1195A Sensor to any finger (except thumb), or to any toe. Finger or toe diameter should measure between 7-8mm (0.27-0.31").

Step	Applying the Glove Style Sensor
1	Select appropriate sensor for patient size (defined above).
2	Position sensor over patient's finger. If applying an adult finger sensor, be sure the sensor cable is positioned ON TOP of the finger/hand as shown above.
3	Patient fingertip should touch but not protrude from end of sensor. If necessary, trim the fingernail to correctly position the sensor.
4	If applying an adult finger sensor, secure the cable to BACK OF HAND with M1627A wristband (supplied only with adult finger sensors.)
5	Plug sensor into instrument (or into adapter cable if needed).
6	Inspect and change the sensor application site periodically.

Applying Wrap Style Sensors

Neonate Patients (1-4kg)

Apply M1193A Sensor to hand or foot as shown.

Step	Applying the Wrap Style Sensor
1	Position the sensor so that it surrounds the hand or foot with optical components opposite each other.
2	Holding sensor, stretch strap slightly (not more than 2.5cm (1.0")).
3	Insert stretched strap into slot, hold it there while threading end through latch. If strap is too long, thread it through second latch.
4	Plug sensor into instrument (or into adapter cable if needed).
5	Inspect and change the sensor application site periodically.

Applying Ear-Clip Style Sensors

Adult or Pediatric Patients (>40kg)

Clip M1194A Sensor to ear as shown.

Step	Applying the Ear-Clip Style Sensor
1	To improve perfusion, massage or warm the ear lobe.
2	Clip probe onto fleshy part of the lobe. Do not position probe on cartilage or pressed against the head. Clip plastic fixing mechanism on top of the ear. The fixing mechanism minimizes artifact generated by patient motion.
3	Plug sensor into instrument (or into adapter cable if needed).
4	Inspect and change the sensor application site periodically.

Cleaning and Low-Level Disinfection

These reusable sensors should be cleaned and disinfected, but never sterilized. Follow the procedure outlined below.

Warnings

- Use only validated cleaning agents and disinfectants listed below; use of others that are not validated may damage the sensor or its connecting wires and shorten product lifetime or cause safety hazards.
- Select disinfectants carefully as some have very similar names but completely different compositions.
- Do not immerse the sensor connector in any of the cleaning solutions, disinfectant, or other liquid (the sensor and cable housing may be immersed, but not the connector).
- Do not soak sensors in disinfectants longer than specified by the disinfectant manufacturer.
- Do not sterilize or autoclave the sensors.

Validated Cleaning Agents

- Mild Detergent
- Salt Solution (1%)

Validated Disinfectants

Metricide [®] 28	Cidex [®] Formula 7	Kohrsolin [®] (2%)
Metricide [®] Plus 30	Cidex [®] OPA	Mucapur [®] -CD (1%)
Terralin [®] Liquid	Cidex [®] Plus	Isopropanol (70%) or Isopropanol Wipe (70%)
Incidin [®] Liquid	Omnicide [®] 28	

Step	Cleaning and Low-Level Disinfection
1	Clean sensor following instructions supplied with cleaning agent.
2	Disinfect sensor following instructions supplied with disinfectant.
3	Rinse sensor and cable (but not the connector) in water, wiping it dry with a clean cloth, leaving to dry completely. If you notice any deterioration or damage to sensor or cable, dispose the sensor immediately.

Specifications

Sensor Accuracy

A-Series Sensors (M1192A, M1193A, M1194A, M1195A): For A-Series sensor SpO₂ and Pulse Rate accuracy specifications, reference your Philips/Agilent/HP instrument's Instructions for Use.

B-Series Sensors (M1191B, M1191BL): SpO₂ and Pulse Rate accuracy specifications for M1191B and M1191BL sensors are the same as indicated for obsolete M1191A and M1191AL sensors in your Philips/Agilent/HP instrument's Instructions for Use. (Note: M1191B and M1191BL sensors replace obsolete M1191A and M1191AL sensors.)

Neonatal Applications: The claimed accuracy numbers for neonatal application are already increased by an additional 1% to account for affect on oximeter measurements due to fetal hemoglobin in neonatal blood, as predicted by the literature.

Light Emitting Diode Wavelength Ranges

Wavelength ranges for the light emitting diodes used in these sensors are within 600nm-1000nm, with an optical output power of less than 15mW. Knowing the wavelength range can be useful to clinicians performing photodynamic therapy.

Measurement Validation

The SpO₂ accuracy has been validated in human studies against arterial blood sample reference measured with a CO-oximeter. In a controlled desaturation study, healthy adult volunteers with saturation levels between 70% and 100% SaO₂ were studied. The population characteristics for those studies were:

- Approximately 50% female and 50% male ranging in age from 18-45
- Skin tone: from light to dark

Because pulse oximeter equipment measurements are statistically distributed, only approximately 2/3 of pulse oximeter equipment measurements can be expected to fall within the $\pm$ Arms value measured by a CO-oximeter. Functional testers, such as an SpO₂ simulator cannot be used to assess the accuracy of pulse oximeter Sensors.

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For more information, please call your local Philips sales office listed in your telephone directory or a Philips regional office listed below for the location of your nearest sales office.

United States:

Philips Medical Systems
North America Corporation
3000 Minuteman Road
Andover, MA 01810-1099
(800) 225-0230

Canada:

Philips Medical Systems Canada
281 Hillmount Road
Markham, ON
L6C 2S3
(800) 291-6743

Latin America Headquarters:

Philips Medical Systems
1550 Sawgrass Corporate Parkway #300
Sunrise, FL 33323 USA
Tel: (954) 835-2600
Fax: (954) 835-2626

Europe, Middle East and Africa:

Philips Medizin Systeme
Böblingen GmbH
Cardiac and Monitoring Systems
Hewlett-Packard Str. 2
71034 Böblingen, Germany
Fax: (+49) 7031 463 1552

Asia Pacific Headquarters:

Philips Medical Systems
30/F Hopewell Centre
17 Kennedy Road
Wanchai
Hong Kong
Tel: (852) 2821 5888
Fax: (852) 2527 6727

 **0366** This product complies with the requirements of the Council Directive 93/42/EEC June 1993 (Medical Device Directive).

<http://www.philips.com/healthcarestore>

Manufacturing Address:

Philips Medizin Systeme
Böblingen GmbH
Hewlett-Packard Str. 2
71034 Böblingen, Germany

Subject to modification
Copyright © 2009
Koninklijke Philips Electronics N.V.
Printed in India March 2009
Philips Part Number 453564154671