




PHILIPS HEALTHCARE	Product Recycling Passport	Page 1 of 2
-------------------------------	-----------------------------------	-------------

Product name:	Certeray iX with front sheets and optional UPS-support	
Identification code(s)	9890-000-65306 and higher 12NC; 9896-010-4030X; optional 4512-105-1617X	
Total weight (in Kg)	184kg and 43,5kg cabinet covers (not pictured); optional 4,5kg for UPS support	
Producer/ Manufacturer	Name company:	Philips Medical Systems DMC GmbH
	Address:	Röntgenstrasse 24
	Zip code:	22335 Hamburg
	Country:	Germany

Recycle Info	Items:	Location
General Hazards (Desinfection)	For dismantling activities Treatment Facilities must consider the national requirements. For personnel which can come into contact with contaminated material, preventive measures pursuant to national requirements must be taken into account. Prior to this instruction, please refer to the System Recycling Passport.	
Special attention 	Items:	
	Some mechanical components are heavier than 25kg	
 Fluids / Gases	Items:	Location
	Cooling Oil, Shell Diala S3 ZX-IG	Figure 1 [1]
	Transformer Oil, Shell Diala	Figure 2 [2]
Batteries  To be Removed	Type:	Location
	N.A.	
Hazardous	Substances:	Location
	All components are RoHS compliant	

 To be Removed		

Locations of materials mentioned on the previous sheet:




Figure 1
Certeray iX front view




Figure 2
Certeray iX side view

Optional:



Figure 3
UPS-support