



System One REMstar SE

An entry-level CPAP device ready for today's challenges

Designed around your needs

As an extension of the System One sleep therapy family, the REMstar SE meets today's CMS reimbursement requirements and provides a robust set of features in one affordable unit. With the REMstar SE, your patients can enjoy many of the same comfort and compliance benefits of System One, including:

- Flex pressure relief
- Enhanced humidification with Heated Tube option
- Downloadable compliance data capable (modem or SD card)
- Enhanced VIC reporting
- EncoreAnywhere compliance management system

Simplicity through operational efficiencies

Philips Respironics provides cost savings and convenience to home care providers by offering one family of devices that share common components such as humidifiers, filters, tubing, and more.

Addresses today's landscape

- Competitive bidding
- Managed care
- Secondary CPAP device



The REMstar SE is just one more way we are working with you as Allies in Better Sleep and Breathing.



Flexing its therapy muscles

Philips Respironics Flex technology has become our clear standard for pressure relief therapy.

A natural way of breathing

Flex technology helps make PAP therapy more natural by its ability to sense when the patient is exhaling and then reducing pressure to make the patient's breathing more natural and comfortable. REMstar SE is covered under our Flex Promise program. If a patient is non-compliant on Flex, we will credit the home care provider 100% of the purchase price toward the purchase of a bi-level OSA device.

Flex



Ordering information

Devices	Part number
REMstar SE	DS220S
REMstar SE with Heated Humidifier	DS220HS
REMstar SE with Heated Tube Humidifier	DS220TS
Accessories	
System One 60 Series Heated Humidifier	DS6H
System One 60 Series Heated Tube Humidifier	DS6T
SD cards (10/pack)	1063859
SD card mailers (10/pack)	1065146
Wireless modem	100200C

Specifications

Pressure range	4 to 20 cm H ₂ O
Flex pressure relief	0 to 3
Ramp time	0 to 45 min (5-min increments)
Starting ramp pressure	4 to CPAP; patient adjustable
Dimensions	7" L x 5.5" W x 4"H 18 cm L x 14 cm W x 10 cm H
Weight	3lbs; 1.36 kg (without humidifier)
Filters	Pollen (foam) and optional ultra-fine
Device setup	LED/keypad
Data storage capacity	Compliance VIC 7-30-day averages SD card > 1 year On-board > 1 year
Compliance meter	Breathing detection
Altitude compensation	Low, medium or high
Electrical requirements	100-240 VAC, 50/60 Hz
Humidification	Integrates with REMstar System One 60 Series Heated Humidifier and System One 60 Series Heated Tube Humidifier
Warranty	Two years
Miscellaneous	Fine pressure adjustment

Respironics, REMstar, System One, and Flex are trademarks of Koninklijke Philips Electronics N.V. and its affiliates. All rights reserved.

Please visit www.philips.com/systemone



© 2013 Koninklijke Philips Electronics N.V. All rights are reserved.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

CAUTION: US federal law restricts these devices to sale by or on the order of a physician.

SB 5/9/13 MCI 4105615 PN 1109871

Philips Healthcare is part of Royal Philips Electronics

Philips Respironics International Headquarters
+33 1 47 28 30 82
Philips Respironics Asia Pacific
+65 6882 5282
Philips Respironics Australia
+61 (2) 9666 4444
Philips Respironics United Kingdom
+44 800 1300 845

Philips Respironics
1010 Murry Ridge Lane
Murrysville, PA 15668

Customer Service
+1 724 387 4000
800 345 6443 (toll free, US only)
www.philips.com/respironics