



# Easy access

**Encore patient management** makes it quick and simple to obtain vital patient therapy intelligence.



Knowing exactly what patients are experiencing at any given moment of therapy is key to making intelligent patient management decisions.

Encore patient management contributes to an intelligent sleep therapy system by providing timely and secure access to comprehensive patient therapy data. Through an intelligent web-based platform, EncoreAnywhere delivers on its name, providing the care team – regardless of where its members are located – with instant access to patient data via new wired or wireless modem options. As a result, the team can collaborate and respond with prescription changes in real time. EncorePro provides the same level of data via an SD card. Both systems are intuitive and secure, making it simpler and more practical to access the patient intelligence that contributes to more intelligent patient and business management.

## **EncoreAnywhere: intelligence at your fingertips**


With the new wireless or wired modem connection, EncoreAnywhere automatically downloads data for instant access. Through a secured-access website available via any computer, physicians, sleep lab professionals, and home care providers can view timely, accurate patient information and then instantly change prescriptions or communicate as necessary. Thanks to a highly intuitive interface, putting such intelligence to work is surprisingly simple, efficient and cost effective. Better patient compliance and increased care team efficiency; that's intelligent therapy management.

## **EncorePro: simple and secure**

The next generation of EncorePro is better than ever for day-to-day patient management. We've simplified the interface while beefing up security, making all of the notifications, reminders, contact records, patient data and reports more secure and even easier to access.

**PHILIPS**

**RESPIRONICS**



An easy-to-navigate **My Day** page helps users prioritize patient intervention needs quickly.

**Prescription** changes can be made securely and as swiftly as the push of a button. Updates to the device are automatic.

Summary and daily detail reports of **Therapy** data – including patient flow waveforms – are easily viewed and arranged by a specific time frame, including our “Best 30” report which provides the best 30 consecutive days of greater than four-hour sessions. In addition, internal device memory stores six months of compliance data and five days of patient flow waveforms that can be transferred by the Encore data card or modem at any time.

Item	Part number
<b>EncoreAnywhere</b>	
EncoreAnywhere subscription	1999999
EncoreAnywhere administrator guide - homecare provider	1056401
EncoreAnywhere administrator guide - sleep lab and physician	1062233
Take Ten referral marketing package	1056397
Take Ten referral set-up and training package	1056398
<b>EncorePro</b>	
EncorePro 2 software	1054785
EncorePro 2 user guide	1060291
EncorePro 2 installation guide	1054258
<b>Accessories</b>	
SD card reader	1047300
SD card cover (attached to device)	1063858
SD cards (10 pack)	1063859
SD card mailer (10 pack)	1065146
Encore reports manual	1055272
<b>Sleep therapy modem systems</b>	
Wired flow modem	100200W
Wired flow modem monthly service	100201W
Wireless flow modem	100200C
Wireless flow modem monthly service	100201C

Respironics, Encore, and EncoreAnywhere are trademarks of Respironics, Inc. and its affiliates.



© 2009 Koninklijke Philips Electronics N.V. All rights are reserved.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

CAUTION: US federal law restricts these devices to sale by or on the order of a physician.

Hoech KW 07/02/09 MCI 4102519 PN 1064785

Philips Healthcare is part of Royal Philips Electronics

Respironics Asia Pacific  
+65 6298 1088  
Respironics Australia  
+61 (2) 9666 4444  
Respironics Europe, Middle East, Africa  
+33 1 47 52 30 00  
Respironics United Kingdom  
+44 800 1300 845

Philips Respironics  
1010 Murry Ridge Lane  
Murrysville, PA 15668

Customer Service  
+1 724 387 4000  
800 345 6443 (toll free, US only)  
www.philips.com/respironics