





## Recycling Passport

<b>Product name:</b>	<b>RO 1750                      SRO 33100                      ROT 380</b>		
<b>Identification code</b>	<b>9890-000-86461   9890-000-86481</b>		
<b>Order No.</b>			
<b>Device Type</b>	<b>B2B (Business to business)</b>		
<b>Producer /Manufacturer</b>	Company name:	Philips Medical Systems DMC GmbH	
	Address:	Roentgenstrasse 24	
	Zip code:	22335 Hamburg	
	Country:	Germany	
	Web info:	<a href="http://www.healthcare.philips.com/main/about/Sustainability/Recycling/">http://www.healthcare.philips.com/main/about/Sustainability/Recycling/</a>	
<b>Total weight (in Kg) thereof</b>			<b>25,627</b>

steel, iron	iron, low alloy (<5%)	1,690	KG
	iron, high alloy (>5%)	1,009	KG
	cast iron, cast steel	0,257	KG
<b>steel, iron</b>		<b>2,956</b>	<b>KG</b>
nonferrous metals and alloys	aluminum, -alloy	4,625	KG
	molybdenum, -alloy	0,567	KG
	copper, -alloy	0,300	KG
	nickel, -alloy	0,134	KG
	tungsten, -alloy	0,088	KG
	zinc, -alloy	0,076	KG
<b>nonferrous metals and alloys</b>		<b>5,790</b>	<b>KG</b>
glass / ceramics	glass	0,341	KG
<b>glass / ceramics</b>		<b>0,341</b>	<b>KG</b>
plastics / organic substances	oil	6,000	KG
	thermoplastic	1,769	KG
	elastomer	0,123	KG
<b>plastics / organic substances</b>		<b>7,892</b>	<b>KG</b>
standard parts	other electronic powered devic	3,040	KG
	mounting parts, attaching part	0,237	KG
	computers and accessories	0,021	KG
<b>standard parts</b>		<b>3,298</b>	<b>KG</b>
relevant materials	lead and -compounds (cables an	5,350	KG
<b>relevant materials</b>		<b>5,350</b>	<b>KG</b>

Printed copies are uncontrolled documents

## Recycling Passport

<b>Recycle info</b>		
<b>General Hazards (Desinfection)</b>	For dismantling activities Treatment Facilities must consider the national requirements. For personnel which can come into contact with contaminated material, preventive measures pursuant to national requirements must be taken into account	
<b>Hazardous Substances</b>    <b>To be removed</b>	<p style="text-align: center;"><b>Substances:</b></p> Lead as x-ray shielding inside housing: 5,35 kg	<p style="text-align: center;"><b>Location</b></p> Figure 1
<b>Batteries</b>    <b>To be removed</b>	<p style="text-align: center;"><b>Type:</b></p>	<p style="text-align: center;"><b>Location</b></p>
<b>Special attention</b>  	<p style="text-align: center;"><b>Item</b></p>	<p style="text-align: center;"><b>Location</b></p>
<b>Fluids / Gases</b>  	<p style="text-align: center;"><b>Item</b></p> Oil: 6 kg  HV-transformer oil contains no PCBs	<p style="text-align: center;"><b>Location</b></p>

Printed copies are uncontrolled documents

# Recycling Passport

Locations as mentioned in the Passport ( Pictures information)



Figure 1

Printed copies are  
uncontrolled documents