




Product name	Set PHD-3/3.1
Identification code(s)	4522 132 78231 / 4522 132 78241
Total weight [kg]	6 kg (incl. packaging)
Manufacturer	Conrac GmbH
Address	Lindenstraße 8
Zip code	D-97990 Weikersheim
Country	Germany
Electronic info	45221327823x_579_v01.doc

Recycling information

Special Attention	Items	Location
		Figure 2, 3

Batteries	Items, type	Location
		

Hazardous substances	Items	Location
	Indicate PCB's here also!	Figure 2

Fluids and gases	Items	Location
		Figure 2

Name: Recycling passport PHD3/3.1-Set		12NC: 4522 132 78231		Version: 00	
Author: Wolfgang Klein	Page 1 of 4	Group No: -	Section: -	Doc.Date: 2008-04-30	
				Status: Released	

Locations as mentioned in the Passport (Information in pictures)

Figure1: Complete PHD3-Set



Materials used:

PHD-Set including TFT-LCD (incl. PCB), Fiber Optic transmitter (incl. PCB), cables
The front bezel of the TFT-LCD is made of Polycarbonate
The housing of the Fiber Optic transmitter is made of aluminium

TFT-LCD (PCB)

Fiber Optic transmitter (PCB)

Cables

Name: Recycling passport PHD3/3.1-Set	12NC: 4522 132 78231	Version: 00 Doc.Date: 2008-04-30
Author: Wolfgang Klein	Page 2 of 4	Group No: - Section: - Status: Released

Figure2 TFT-LCD



Mercury contained in
backlight hidden behind
LCD screen

Figure3 Fiber Optic transmitter




Name: Recycling passport PHD3/3.1-Set		12NC: 4522 132 78231		Version: 00 Doc.Date: 2008-04-30	
Author: Wolfgang Klein		Page 3 of 4	Group No: -	Section: -	Status: Released

Recycling Passport Revision History

Revision	Description	Date
1	Initial version	25-04-2008

Recycling Passport Approval Signatures

<p>Written By</p>  <p>CONRAC GmbH (Wolfgang Klein)</p>	<p>Reviewed By (optional)</p>
<p>Approved By</p> <p>D.Sinnema</p>	<p>Approved By</p>

End of document.

Name: Recycling passport PHD3/3.1-Set		12NC: 4522 132 78231	Version: 00 Doc.Date: 2008-04-30
Author: Wolfgang Klein	Page 4 of 4	Group No: -	Section: - Status: Released