



<b>Model 274 Gradient Amplifier System</b>	<b>Copley part number 00-00272-000</b>
Identification code(s)	4522-131-54201
Total weight [kg]	696.6 kg
<b>Manufacturer</b>	<b>Copley Controls Corporation</b>
Address	20 Dan Road, Canton Massachusetts
Zip code	02021
Country	USA
Electronic info	Gradient Amplifier C274_4522-131-5420x.pdf

## Recycling information

<b>Special Attention</b>	<b>Items</b>	<b>Location</b>
	PVC wiring & cabling	Internal to Rack
	Chromate finish (all chassis)	Internal to Rack
	Nickle plated hardware	Internal to Rack

<b>Hazardous substances</b>	<b>Items</b>	<b>Location</b>
	Printed circuit boards utilize lead (all assemblies)	Internal to Rack

<b>Fluids and gases</b>	<b>Items</b>	<b>Location</b>
		

Name: <b>Recycling passport Copley 274</b>		12NC: <b>4522 131 54200</b>		Version: <b>02</b> Doc.Date:09-Aug-2007	
Author:		Pag.: <b>1 of 1</b>	Group No: 5790	Section: -	Status: <b>Authorized</b>