






<b>Product name:</b>	OP.CONS. +OCD-XW8200-C-M		
<b>Identification code(s)</b>	9896 030 14231		
<b>Total weight (in Kg)</b>	53 Kg Excluding Computer , Display and DVD writer		
<b>Producer /Manufacturer</b>	Name company:	Philips Medical Systems Nederland B.V.	
	Address:	Veenpluis 4-6	
	Zip code:	5680 DA BEST	
	Country:	The Netherlands	
	Electronic info:		

<b>Recycle info</b>	<b>Substances:</b>	<b>Location</b>
<b>Special attention</b> 	<b>Items:</b>	<b>Location</b>
	Contains PWB materials.	Figure 1 and 5
	Computer XW8200 not part of this passport see 4522 132 44812	Figure 2
	Display OCD-set not part of this passport see 4522 132 34592	Figure 3
	DVD writer not part of this passport see 4522 132 33791	Figure 4
<b>Batteries</b>  <b>To be removed</b>	<b>Type:</b>	<b>Location</b>
	NA	
<b>Hazardous</b>  <b>To be removed</b>	<b>Examples</b>	
	NA	
 <b>Fluids / Gases</b>	<b>Items:</b>	<b>Location</b>
	NA	

Name: <b>Recycling Passport Console</b>		12NC: <b>9896 030 14230</b>		Version: <b>01</b> Doc.Date: <b>2007-07-31</b>	
Author:		Pag.: <b>1 of 2</b>	Group No: <b>5790</b>	Section: <b>-</b>	Status: <b>Authorized</b>



**Figure 2**




**Figure 1**



Keyboard, mouse



**Figure 3**



Display



**Figure 4**



**Figure 5**

Name: <b>Recycling Passport Console</b>		12NC: <b>9896 030 14230</b>		Version: <b>01</b> Doc.Date: <b>2007-07-31</b>	
Author:		Pag.: <b>2 of 2</b>	Group No: <b>5790</b>	Section: <b>-</b>	Status: <b>Authorized</b>