





Product name:	dS Head coil 3.0T		
Identification code(s)	4598 000 9798x		
Total weight (in Kg)	3.0 Kg		
Producer /Manufacturer	Name company:	Philips Healthcare Best	
	Address:	Veenpluis 4-6	
	Zip code:	5680DA Best	
	Country:	The Netherlands	
	Electronic info:	NVC Head 3.0T	

Recycle info	Substances:	Location
Special attention 	Items:	
	26 pcb assemblies and one flexfoil pcb assembly	Figure 1
	2 so-called RXE modules: 45221330827x or 45980004622x (containing pcb assy)	Figure 2
 Batteries <i>To be removed</i>	Type:	Location
	No batteries present.	
Hazardous  To be removed		Location
	SMP connectors (80 total) contain BeCu	Figure 3
 Fluids / Gases	Items:	Location
	No fluids and/or gases present.	

Name:	WEEE NVC Head 3.0T	12NC:	DMR152837	Version: 00	Doc.Date: 2012-07-31
Author:		Pag.: 1 of 4	Group No: 5790	Section: -	Status: see Agile

Locations as mentioned in the Passport (Picture information)



Figure 1, NVC Head 3.0T pcb locations

Name: WEEE NVC Head 3.0T		12NC: DMR152837		Version: 00 Doc.Date: 2012-07-31	
Author:		Pag.: 2 of 4	Group No: 5790	Section: -	Status: see Agile




Figure 2, NVC Head 3.0T RXE locations

Name: WEEE NVC Head 3.0T		12NC: DMR152837		Version: 00 Doc.Date: 2012-07-31	
Author:		Pag.: 3 of 4	Group No: 5790	Section: -	Status: see Agile




Figure 3, SMP Connector (multiple locations) in NVC Head 3.0T

Name: WEEE NVC Head 3.0T		12NC: DMR152837		Version: 00 Doc.Date: 2012-07-31	
Author:		Pag.: 4 of 4	Group No: 5790	Section: -	Status: see Agile
© PHILIPS MEDICAL SYSTEMS NEDERLAND BV, BEST 2012				Template version 2.	