




Product name:	PFEI-QBC-iRF 15T/30T		
Identification code(s)	9896-030-2569 x / 9896-030-2602 x		
Total weight (in Kg)	20 Kg		
Producer /Manufacturer	Name company:	Philips Healthcare Best	
	Address:	Veenpluis 4-6	
	Zip code:	5680DA Best	
	Country:	The Netherlands	
	Electronic info:	PFEI-QBC-iRF 15T/30T	

Recycle info	Substances:	Location
Special attention 	Items:	
	No special attention needed.	
 Batteries To be removed	Type:	Location
	No batteries present.	
Hazardous  To be removed	Examples	Location
	Depending on the implementation four or six Printed circuit board assemblies	Figure 1a and 1b
	Each Printed circuit board assembly contains 10 times 14mm beryllium copper sticky fingers.	Figure 2
 Fluids / Gases	Items:	Location
	No fluids and/or gases present.	

Name: WEEE PFEI-QBC-iRF 15T/30T, 9896-030-2569 / 9896-030-2602		12NC: DMR162786		Version: 00.1 Doc.Date: 2013-01-14	
Author:		Pag.: 1 of 2	Group No: 5790	Section: -	Status: Authorized

Locations as mentioned in the Passport (Pictures information)




Figure 1a, PFEI-WA-iRF right-side




Figure 1b, PFEI-WA-iRF left-side




Figure 2, Beryllium Sticky-Fingers on Printed circuit board assembly

Note: Figures comprehends the none MN option

Name: WEEE PFEI-QBC-iRF 15T/30T, 9896-030-2569 / 9896-030-2602		12NC: DMR162786		Version: 00.1 Doc.Date: 2013-01-14	
Author:		Pag.: 2 of 2	Group No: 5790	Section: -	Status: Authorized