





Product name	USER INTERFACE MODULE RIGHT & USER INTERFACE MODULE LEFT TRX
Identification code(s)	4522 132 67641 & 4522 132 76241
Total weight [kg]	1,780 kg (1,052 kg - UI only) & 1,781 kg (1,053 kg – UI only)
Manufacturer	Philips Medical Systems Nederland B.V.
Address	Veenpluis 4-6, Best
Zip code	5684 PC
Country	The Netherlands
Electronic info	

Recycling information

Special Attention	Items	Location
	Modules contains a substantial amount of V2A Steel (X5 CRNI 18-10)	Figure1

Batteries	Items, type	Location
	No batteries contained in the product.	

Hazardous substances	Items	Location
	Printed Circuit Board is manufactured according RoHS directive.	Figure 2

Fluids and gases	Items	Location
	No fluids or gases contained in the product.	

Name: Recycling Passport USER INTERFACE MODULE RIGHT		12NC: 4522 132 67640		Version: 00 Doc.Date: 2009-08-31	
Author:		Pag.: 1 of 2	Group No: 5790	Section: -	Status: Approved

Locations as mentioned in the Passport (Information in pictures)



Figure1





Figure2



Name: Recycling Passport USER INTERFACE MODULE RIGHT		12NC: 4522 132 67640		Version: 00 Doc.Date: 2009-08-31	
Author:		Pag.: 2 of 2	Group No: 5790	Section: -	Status: Approved