

Product name:	Anterior Coil M2 1.5T		
Identification code(s)	4598-001-1839x		
Total weight (in Kg)	5.1 Kg		
Producer /Manufacturer	Name company:	Philips Healthcare Best	
	Address:	Veenpluis 4-6	
	Zip code:	5680DA Best	
	Country:	The Netherlands	
	Electronic info:	Anterior Coil M2 1.5T	

Recycle info	Substances:	Location
Special attention 	Items:	
	Sixteen PCB assemblies and one Flex-foil PCB assembly	Figure 1
	Four so-called RXE modules: 4598-000-4622x (containing PCB assy)	Figure 2
 Batteries To be removed	Type:	Location
	No batteries present.	
Hazardous  To be removed		Location
	SMP connectors (80 total) contain BeCu	Figure 3
 Fluids / Gases	Items:	Location
	No fluids and/or gases present.	

Name: WEEE Anterior Coil M2 1.5T		12NC: DMR149398		Version: 00 Doc.Date: 2012-06-11	
Author:		Pag.: 1 of 4	Group No: 5790	Section: -	Status: see Agile

Locations as mentioned in the Passport (Picture information)




Figure 1. Anterior Coil M2 1.5T pcb locations

Name: WEEE Anterior Coil M2 1.5T		12NC: DMR149398		Version: 00 Doc.Date: 2012-06-11	
Author:		Pag.: 2 of 4	Group No: 5790	Section: -	Status: see Agile
© PHILIPS MEDICAL SYSTEMS NEDERLAND BV, BEST 2012				Template version 2.	




Figure 2. Anterior Coil M2 1.5T RXE locations

Name: WEEE Anterior Coil M2 1.5T		12NC: DMR149398		Version: 00 Doc.Date: 2012-06-11	
Author:		Pag.: 3 of 4	Group No: 5790	Section: -	Status: see Agile
© PHILIPS MEDICAL SYSTEMS NEDERLAND BV, BEST 2012				Template version 2.	




Figure 3. SMP Connector (multiple locations) in AC M2 Assembly

Name: WEEE Anterior Coil M2 1.5T		12NC: DMR149398		Version: 00 Doc.Date: 2012-06-11	
Author:		Pag.: 4 of 4	Group No: 5790	Section: -	Status: see Agile
© PHILIPS MEDICAL SYSTEMS NEDERLAND BV, BEST 2012				Template version 2.	