

Product name:	NVC Head Neck 1.5T		
Identification code(s)	4598-001-7123x		
Total weight (in Kg)	3.4 Kg		
Producer /Manufacturer	Name company:	Philips Healthcare Best	
	Address:	Veenpluis 4-6	
	Zip code:	5680DA Best	
	Country:	The Netherlands	
	Electronic info:	NVC Head Neck 1.5T	

Recycle info	Substances:	Location
Special attention 	Items:	
	No special attention needed.	
 Batteries To be removed	Type:	Location
	No batteries present.	
Hazardous  To be removed		Location
	Eight pcb assemblies and two flexfoil pcb assemblies	Figure 1
	Two so-called RXE modules: 45221330827x or 45980004622x (containing pcb assy)	Figure 2a, 2b
	Twenty-four copper-gold contact pins containing BeCu springs	Figure 3
 Fluids / Gases	Items:	Location
	No fluids and/or gases present.	

Name: WEEE NVC Headneck 1.5T	12NC: DMR148346	Version: 01 Doc.Date: 2012-05-22
Author:	Pag.: 1 of 3	Group No: 5790
	Section: -	Status: see Agile

Locations as mentioned in the Passport (Picture information)




Figure 1, NVC Head Neck 1.5T pcb locations




Figure 2a, 2b, NVC Head Neck 1.5T RXE locations

Name: WEEE NVC Headneck 1.5T	12NC: DMR148346	Version: 01 Doc.Date: 2012-05-22
Author:	Pag.: 2 of 3	Group No: 5790
		Section: -
		Status: see Agile
© PHILIPS MEDICAL SYSTEMS NEDERLAND BV, BEST 2011		Template version 2.




Figure 3, NVC Head Neck 1.5T contact pin locations (protruding bottom cover)

Name: WEEE NVC Headneck 1.5T		12NC: DMR148346		Version: 01 Doc.Date: 2012-05-22	
Author:		Pag.: 3 of 3	Group No: 5790	Section: -	Status: see Agile
© PHILIPS MEDICAL SYSTEMS NEDERLAND BV, BEST 2011				Template version 2.	