





<b>Product name</b>	<b>OCD4 (OPERATOR CONSOLE MONITOR-FIMI-Q19)</b>
Identification code(s)	4522 132 3316x which is part of OCD-set (incl.cables) 4522 132 34592
Total weight [kg]	8 kg
<b>Manufacturer</b>	<b>FIMI</b>
Address	Via Saul Banfi, 1 - 21047 Saronno
Zip code	
Country	Italy
Electronic info	

## Recycling information

Special Attention	Items	Location
	Printed Circuit Boards (produced by lead free soldering process)	Figure 1 (item2)

Batteries	Items, type	Location
	NO	Figure ..

<b>Recycle info</b> <b>Hazardous</b>  <b>To be removed</b>	Substances		Location
	Type	Quantity	
	Cd	0	
	Hg	2,5mg max. per lamp(*)	Figure 1 ( item 1 )
	Pb	0	
	Cr <sup>6+</sup>	0	
	PBB	0	
	PBDE	0	
	(*) Mercury is present in the 4 lamps of Backlight, according to RoHS limitation		

Fluids and gases	Items	Location
	NO	Figure ..

Name:	<b>Recycling passport OCD 4 - FIMI-Q19</b>	12NC:	<b>4522 132 33160</b>	Version: <b>02</b>	Doc.Date: <b>2008-09-25</b>
Author:		Pag.: <b>1 of 2</b>	Group No: <b>5790</b>	Section: <b>-</b>	Status: <b>Approved</b>

Locations as mentioned in the Passport ( Pictures information)




FIG. 1

Monitor material		
Fe	620 g	
Al	590 g	
Cu	100 g	
Plastics	1100 g	
Boards (S <sup>2</sup> > 10cm <sup>2</sup> )	98 cm <sup>2</sup> / 65 g	Scaler Board
	35 cm <sup>2</sup> / 15 g	Keyboard
	225 cm <sup>2</sup> / 270 g	Power Board
LCD Panel	2980 g	19"
Pedestal material		
Fe	1920 g	
Plastics	480 g	

Name:	<b>Recycling passport OCD 4 - FIMI-Q19</b>	12NC:	<b>4522 132 33160</b>	Version: <b>02</b>
Author:		Group No: 5790	Section: -	Doc.Date: <b>2008-09-25</b>
	Pag.: <b>2 of 2</b>			Status: <b>Approved</b>