





Product name	OCD4 (OPERATOR CONSOLE MONITOR-FIMI-Q19)
Identification code(s)	4522 132 3316x which is part of OCD-set (incl.cables) 4522 132 34592
Total weight [kg]	8 kg
Manufacturer	FIMI
Address	Via Saul Banfi, 1 - 21047 Saronno
Zip code	
Country	Italy
Electronic info	

Recycling information

Special Attention 	Items	Location
	Printed Circuit Boards (produced by lead free soldering process)	Figure 1 (item2)


Batteries 	Items, type	Location
	NO	Figure ..

Recycle info Hazardous  To be removed	Substances		Location
	Type	Quantity	
	Cd	0	
	Hg	2,5mg max. per lamp(*)	Figure 1 (item 1)
	Pb	0	
	Cr ⁶⁺	0	
	PBB	0	
	PBDE	0	
(*) Mercury is present in the 4 lamps of Backlight, according to RoHS limitation			

Fluids and gases 	Items	Location
	NO	Figure ..

Name:	Recycling passport OCD 4 - FIMI-Q19	12NC:	4522 132 33160	Version: 02
Author:		Group No: 5790	Section: -	Doc.Date: 2008-09-25
	Pag.: 1 of 2			Status: Approved

Locations as mentioned in the Passport (Pictures information)



Q19-S series
Q19-L series

FIG. 1

Monitor material		
Fe	620 g	
Al	590 g	
Cu	100 g	
Plastics	1100 g	
Boards (S ² > 10cm ²)	98 cm ² / 65 g	Scaler Board
	35 cm ² / 15 g	Keyboard
	225 cm ² / 270 g	Power Board
LCD Panel	2980 g	19"
Pedestal material		
Fe	1920 g	
Plastics	480 g	

Name: Recycling passport OCD 4 - FIMI-Q19		12NC: 4522 132 33160	Version: 02 Doc.Date: 2008-09-25
Author:	Pag.: 2 of 2	Group No: 5790	Section: - Status: Approved