

| | | | |
|-------------------------------|------------------|-------------------------|--|
| Product name: | NVC Base 3T | | |
| Identification code(s) | 4598-001-7117x | | |
| Total weight (in Kg) | 4.2 Kg | | |
| Producer /Manufacturer | Name company: | Philips Healthcare Best | |
| | Address: | Veenpluis 4-6 | |
| | Zip code: | 5680DA Best | |
| | Country: | The Netherlands | |
| | Electronic info: | NVC Base 1.5T | |

| Recycle info | Substances: | Location |
|---|---|-----------------|
| Special attention  | Items: | |
| | No special attention needed. | |
| | | |
| | | |
|  Batteries To be removed | Type: | Location |
| | No batteries present. | |
| | | |
| | | |
| | | |
| Hazardous  To be removed | | Location |
| | Eleven pcb assemblies and one flexfoil pcb assembly | Figure 1 |
| | Two so-called RXE modules: 45221330827x or 45980004622x (containing pcb assy) | Figure 2 |
| | Eight so-called RF traps (containing flexfoil pcb material) | Figure 3 |
| | Twelve copper-beryllium spring contacts (mounted under protective cover) | Figure 4 |
|  Fluids / Gases | Items: | Location |
| | No fluids and/or gases present. | |
| | | |
| | | |

| | | | | |
|---------|------------------|----------------|-------------------|----------------------|
| Name: | WEEE NVC Base 3T | 12NC: | DMR134426 | Version: 00 |
| Author: | | Group No: 5790 | Section: - | Doc.Date: 2011-10-07 |
| | | Pag.: 1 of 3 | Status: see Agile | |

Locations as mentioned in the Passport (Picture information)




Figure 1, NVC Base 3T pcb locations




Figure 2, NVC Base 3T RXE locations

| | | | | |
|---|--------------------|-------------------------------------|------------|---------------------|
| Name: WEEE NVC Base 3T | 12NC: DMR134426 | Version: 00 Doc.Date: 2011-10-07 | | |
| Author: | Pag.: 2 of 3 | Group No: 5790 | Section: - | Status: see Agile |
| © PHILIPS MEDICAL SYSTEMS NEDERLAND BV, BEST 2011 | | | | Template version 2. |




Figure 3, NVC Base 3T, examples of RF trap locations




Figure 4, NVC Base 3T CuBe spring contact locations (under hatches)

| | | | | | |
|---|--|--------------------|----------------|-------------------------------------|-------------------|
| Name: WEEE NVC Base 3T | | 12NC: DMR134426 | | Version: 00 Doc.Date: 2011-10-07 | |
| Author: | | Pag.: 3 of 3 | Group No: 5790 | Section: - | Status: see Agile |
| © PHILIPS MEDICAL SYSTEMS NEDERLAND BV, BEST 2011 | | | | Template version 2. | |