





Product name:	PFEI-DDAS 15T 8CH	9896-030-2724.x
	PFEI-DDAS 15T 16/32CH	9896-030-2730.x
	PFEI-DDAS 30T 16/32CH	9896-030-2727.x
Total weight (in Kg)	20 Kg	

Producer /Manufacturer	Name company:	Philips Healthcare Best
	Address:	Veenpluis 4-6
	Zip code:	5680DA Best
	Country:	The Netherlands
	Electronic info:	PFEI-DDAS 15T/30T – 8ch/16ch

RECYCLE INFO	DESCRIPTION	LOCATION
<p><i>Special attention</i></p> 	No special attention needed.	n.a.
<p><i>Batteries</i></p>  <p><i>To be removed</i></p>	No batteries present.	n.a.
<p><i>Hazardous</i></p>  <p>To be removed</p>	Depending on the PFEI-DDAS version five or six Printed circuit board assemblies are present. Each Printed circuit board assembly contains 10 times 14mm beryllium copper sticky fingers.	Figure 1
<p><i>Fluids / Gases</i></p> 	No fluids and/or gases present.	n.a.

Name: WEEE PFEI-DDAS	12NC: DMR195004	Version: 00.1 Doc.Date: 15-01-2014
Author:	Pag.: 1 of 2	Group No: 5790 Section: - Status: Authorized

Locations as mentioned in the Passport (Pictures information)

PCB-Assemblies

PCB-Assemblies

*Berillium-Copper
Sticky-Fingers*




Figure 1, PCB-Assemblies in the front and back of the PFEI-DDAS unit and the position of the Sticky-Fingers on the PCB-Assembly

Name: WEEE PFEI-DDAS		12NC: DMR195004		Version: 00.1 Doc.Date: 15-01-2014	
Author:		Pag.: 2 of 2	Group No: 5790	Section: -	Status: Authorized