




Product name:	Physiotrak Left docking Assy Physiotrak Right docking Assy	
Identification code(s)	4522 132 14280 4522 132 14290	
Total weight (in Kg)	200	
Producer /Manufacturer	Name company:	Philips Medical Systems Nederland B.V.
	Address:	Veenpluis 4-6
	Zip code:	5680 DA BEST
	Country:	The Netherlands
	Electronic info:	

Recycle info	Substances:	Location
Special attention 	Items:	Location
	MEB contains 7 Printed Circuit Board's See doc 452213223380 5790	See fig 1
	MDCU contains Printed Circuit Board's See doc 452213223390 5790	See fig 1

Batteries	Type:	Location
 To be removed	2 x Pb battery	See figure 2

Hazardous		
 To be removed		

	Items:	Location
 Fluids / Gases	Not Applicable	

Name:	WEEE Passport Physiotrak	12NC:	4522 132 14280	Version:	00
Author:		Group No:	5790	Doc.Date:	2006-09-04
		Pag.:	1 of 2	Section:	-
				Status:	Authorized



Figure 1



Figure 2

Name: WEEE Passport Physiotrak		12NC: 4522 132 14280		Version: 00 Doc.Date: 2006-09-04	
Author:		Pag.: 2 of 2	Group No: 5790	Section: -	Status: Authorized