


<b>Product name:</b>	Anterior Coil 3T		
<b>Identification code(s)</b>	4510-010-1290x		
<b>Total weight (in Kg)</b>	4.5 Kg		
<b>Producer /Manufacturer</b>	Name company:	Philips Healthcare Best	
	Address:	Veenpluis 4-6	
	Zip code:	5680DA Best	
	Country:	The Netherlands	
	Electronic info:	Anterior Coil 3T	

<b>Recycle info</b>	<b>Substances:</b>	<b>Location</b>
<b>Special attention</b> 	<b>Items:</b>	
	No special attention needed.	
 <b>Batteries</b> <b>To be removed</b>	<b>Type:</b>	<b>Location</b>
	No batteries present.	
<b>Hazardous</b>  <b>To be removed</b>		<b>Location</b>
	Sixteen pcb assemblies and one Flex-foil pcb assembly	Figure 1
	Four so-called RXE modules: 45221330828x or 45980004624x (containing pcb assy)	Figure 2
 <b>Fluids / Gases</b>	<b>Items:</b>	<b>Location</b>
	No fluids and/or gases present.	


Name:	WEEE Anterior Coil 3T	12NC:	DMR125942	Version: 00
Author:		Group No: 5790	Section: -	Doc.Date: 2010-11-16
	Pag.: 1 of 3			Status: see Agile

**Locations as mentioned in the Passport (Picture information)**



*Figure 1, Anterior Coil 3T pcb locations*

Name: WEEE Anterior Coil 3T		12NC: DMR125942		Version: 00 Doc.Date: 2010-11-16	
Author:		Pag.: 2 of 3	Group No: 5790	Section: -	Status: see Agile
© PHILIPS MEDICAL SYSTEMS NEDERLAND BV, BEST 2011				Template version 2.	



*Figure 2, Anterior Coil 3T RXE locations*

Name: WEEE Anterior Coil 3T		12NC: DMR125942		Version: 00 Doc.Date: 2010-11-16	
Author:		Pag.: 3 of 3	Group No: 5790	Section: -	Status: see Agile
© PHILIPS MEDICAL SYSTEMS NEDERLAND BV, BEST 2011				Template version 2.	