





<b>Product name:</b>	dS Base coil 1.5T		
<b>Identification code(s)</b>	4598 000 9794x		
<b>Total weight (in Kg)</b>	4.7 Kg		
<b>Producer /Manufacturer</b>	Name company:	Philips Healthcare Best	
	Address:	Veenpluis 4-6	
	Zip code:	5680DA Best	
	Country:	The Netherlands	
	Electronic info:	NVC Base 1.5T	

<b>Recycle info</b>	<b>Substances:</b>	<b>Location</b>
<b>Special attention</b> 	<b>Items:</b>	
	26 pcb assemblies and one flexfoil pcb assembly	Figure 1
	2 so-called RXE modules: 45221330827x or 45980004622x (containing pcb assy)	Figure 2
 <b>Batteries</b> <b>To be removed</b>	<b>Type:</b>	<b>Location</b>
	No batteries present.	
<b>Hazardous</b>  <b>To be removed</b>		<b>Location</b>
	SMP connectors (80 total) contain BeCu	Figure 3
 <b>Fluids / Gases</b>	<b>Items:</b>	<b>Location</b>
	No fluids and/or gases present.	

Name:	WEEE NVC Base 1.5T	12NC:	DMR152797	Version: 00	Doc.Date: 2012-07-31
Author:		Pag.: 1 of 4	Group No: 5790	Section: -	Status: see Agile

Locations as mentioned in the Passport (Picture information)






Figure 1, NVC Base 1.5T pcb locations

Name: WEEE NVC Base 1.5T		12NC: DMR152797		Version: 00 Doc.Date: 2012-07-31	
Author:		Pag.: 2 of 4	Group No: 5790	Section: -	Status: see Agile
© PHILIPS MEDICAL SYSTEMS NEDERLAND BV, BEST 2012				Template version 2.	



*Figure 2, NVC Base 1.5T RXE locations*

Name: <b>WEEE NVC Base 1.5T</b>		12NC: <b>DMR152797</b>		Version: 00 Doc.Date: 2012-07-31	
Author:		Pag.: <b>3 of 4</b>	Group No: 5790	Section: -	Status: see Agile
© PHILIPS MEDICAL SYSTEMS NEDERLAND BV, BEST 2012				Template version 2.	



*Figure 3, SMP Connector (multiple locations) in NVC Base 1.5T*

Name: WEEE NVC Base 1.5T		12NC: DMR152797		Version: 00 Doc.Date: 2012-07-31	
Author:		Pag.: 4 of 4	Group No: 5790	Section: -	Status: see Agile
© PHILIPS MEDICAL SYSTEMS NEDERLAND BV, BEST 2012				Template version 2.	