





<b>Product name:</b>	DuraDiagnost	
<b>Identification code(s)</b>	989701002021	
<b>Total weight (in Kg)</b>	927	
<b>Producer/ Manufacturer</b>	<b>Name company:</b>	Philips Healthcare (Suzhou) Co., Ltd.
	<b>Address:</b>	No. 258, ZhongYuan Road, Suzhou Industrial Park, Suzhou, Jiangsu Province
	<b>Zip code:</b>	215024
	<b>Country:</b>	CHINA

Recycle Info	Items:	Location
<b>Special attention</b>  		

Fluids / Gases	Items:	Location
	Oil tank / Remove oil before dismantling	1.3(Tube)

Batteries	Type:	Location
  <b>To be Removed</b>	Batteries, Lithium	2(Workstation Tower)
		2 (Wireless detector)
		1 (Table)

Hazardous	Substances:	Location
  <b>To be Removed</b>	Lead	1.1(Detector); 1.2(Tube & Collimator)
	Thallium	1.1(Detector)

**Materials Report**

Geometry & Image Chain Components			
Class	Feature Description	Gross Weight	Unit
Steel, Iron	Iron, low alloy (<5%)	704.682	Kg
	Iron, High alloy (>5%)	32.948	Kg
		737.630	Kg
Nonferrous metals and alloys	Aluminum, -alloy	98.438	Kg
	Copper, -alloy	2.371	Kg
	Zinc, -alloy	0.001	Kg
	Molybdenum, -alloy	0.500	Kg
	111.310	Kg	
Glass/ceramics	Glass	4.250	Kg
		4.250	Kg
Plastics/organic substances	Thermoplastic	9.273	Kg
	Oil	5.802	Kg
	Thermoset	4.187	Kg
	Elastomer	8.895	Kg
	28.151	Kg	
Standard parts	Bearing plates	43	Kg
	cables	77.461	Kg
	Mounting parts, attaching part	1.538	Kg
	Power supplies/transformers	6.268	Kg
	Other electronic powered device	0.998	Kg
	Printed circuit boards	9.459	Kg
	Rasters	5.098	Kg
	Wheels/castors	1.300	Kg
	Motor, pumps	10.2	Kg
	Brake/clutches	0.83	Kg
	156.722	Kg	
Relevant materials	Lead and –compounds	19.868	Kg
	Cobalt and –compounds	0.240	Kg
	Thallium and - compounds	0.001	Kg
	Tungsten and –compounds	0.080	Kg
	Confirmed by supplier	0.000	Kg
	22.389	Kg	
Prohibited materials	Confirmed by supplier	0.000	Kg
		0.000	Kg
		832.863	Kg


Eleva workspot PC			
Class	Feature Description	Gross weight	Unit
Steel, iron	Iron, low alloy (<5%)	3.020	Kg
		3.020	Kg
Batteries	Batteries, lithium	0.003	Kg
		0.003	Kg
Standard parts	Monitors	16.800	Kg
	Computers and accessories	15.975	Kg
	Printed circuit boards	1.230	Kg
	Cables	1.200	Kg
	Other electronic powered devices	1.185	Kg
	Power supplies/transformers	0.865	Kg
	37.255	Kg	
Relevant materials	Confirmed by suppliers	0.000	Kg
		0.000	Kg
Prohibited materials	Confirmed by suppliers	0.000	Kg

Title: \_\_\_\_\_  
 Recycling Passport Number: \_\_\_\_\_ Rev: \_\_\_\_\_


Copies are uncontrolled

	0.000	Kg
	40.278	Kg

Without M-Cabinet, see M-Cabinet Recycling Passport.



**Figure 1**



**Figure 2**

Title: \_\_\_\_\_  
 Recycling Passport Number: \_\_\_\_\_ Rev: \_\_\_\_\_

**Copies are uncontrolled**