


Product name	gMDU
Identification code(s)	9896 030 1995 x 9896 030 2014 x 9896 030 2048.x
Total weight [kg]	90 kg (9896 030 1995 x) 125 kg (9896 030 2014 x) 275 kg (9896 030 2048.x)
Manufacturer	b-NRGY
Address	Collse Heide 48 Nuenen
Zip code	5674VN
Country	The Netherlands
Electronic info	

Recycling information

Special Attention	Items	Location
	19 kVA auto transformer; 35 kg	See figure 1
	30 kVA isolation transformer; 155 kg	See figure 2
	(heavy transformers)	

Batteries	Items, type	Location
	None	

Hazardous substances	Items	Location
	None	

Fluids and gases	Items	Location
	None	

Name: Recycling passport gMDU		12NC: 9896 030 1995.x	Version: 00 Doc.Date: 2010-11-20	
Author:	Page 1 of 3	Group No: 5790	Section: -	Status: Auth.

Locations as mentioned in the Passport (Information in pictures)

Figure 1 gMDU 400V-60Hz; 9896 030 2014.x



The cabinet of the 9896 030 1995.x (gMDU 400V-50Hz) and 9896 030 2014.x (gMDU 480V 60 Hz) are similar except for the 19 kVA transformer.

Figure 2: gMDU 480V-60Hz ; 9896 030 2048.x

See next page

Name: Recycling passport gMDU	12NC: 9896 030 1995.x	Version: 00 Doc.Date: 2010-11-20
Author:	Page 2 of 3	Group No: 5790 Section: - Status: Auth.
© PHILIPS MEDICAL SYSTEMS NEDERLAND BV, BEST 2011		



Name: Recycling passport gMDU		12NC: 9896 030 1995.x		Version: 00 Doc.Date: 2010-11-20	
Author:		Page 3 of 3	Group No. 5790	Section: -	Status: Auth.