


<b>Product name</b>	<b>gMDU</b>
Identification code(s)	9896 030 1995 x 9896 030 2014 x 9896 030 2048.x
Total weight [kg]	90 kg (9896 030 1995 x)  125 kg (9896 030 2014 x)  275 kg (9896 030 2048.x)
<b>Manufacturer</b>	<b>b-NRGY</b>
Address	Collse Heide 48 Nuenen
Zip code	5674VN
Country	The Netherlands
Electronic info	

### Recycling information

Special Attention	Items	Location
	19 kVA auto transformer; 35 kg	See figure 1
	30 kVA isolation transformer; 155 kg	See figure 2
	(heavy transformers)	

Batteries	Items, type	Location
	None	

Hazardous substances	Items	Location
	None	

Fluids and gases	Items	Location
	None	

Name: <b>Recycling passport gMDU</b>		12NC: <b>9896 030 1995.x</b>	Version: <b>00</b> Doc.Date: <b>2010-11-20</b>
Author:	Page <b>1</b> of <b>3</b>	Group No: <b>5790</b>	Section: <b>-</b> Status: <b>Auth.</b>

## Locations as mentioned in the Passport (Information in pictures)

Figure 1 gMDU 400V-60Hz; 9896 030 2014.x



The cabinet of the 9896 030 1995.x (gMDU 400V-50Hz) and 9896 030 2014.x (gMDU 480V 60 Hz) are similar except for the 19 kVA transformer.

Figure 2: gMDU 480V-60Hz ; 9896 030 2048.x

See next page

Name: <b>Recycling passport gMDU</b>		12NC: <b>9896 030 1995.x</b>		Version: <b>00</b> Doc.Date: <b>2010-11-20</b>	
Author:		Page <b>2</b> of <b>3</b>	Group No. <b>5790</b>	Section: <b>-</b>	Status: <b>Auth.</b>
© PHILIPS MEDICAL SYSTEMS NEDERLAND BV, BEST 2011					



Name: <b>Recycling passport gMDU</b>		12NC: <b>9896 030 1995.x</b>		Version: <b>00</b> Doc.Date: <b>2010-11-20</b>	
Author:		Page <b>3 of 3</b>	Group No. <b>5790</b>	Section: <b>-</b>	Status: <b>Auth.</b>