




Product name	DACC(2)
Identification code(s)	See list in below tabel / 4522 133 0896x(Basic cabinet)
Total weight [kg]	260 ± 20 kg
Manufacturer	Philips Medical Systems Nederland B.V.
Address	<i>Veenpluis 4-6 Best</i>
Zip code	5680DA
Country	The Netherlands
Electronic info	DMR192492 WEEE DACC(2)

Recycling information

Special Attention 	1/ DACC2 Basic cabinet 4522 133 0896x	Special attention - NA
	Items: The following items have their own passport and are not part of this passport	Location
	2/ AIBo 4598 001 3851x	Figure 1 DACC(2)
	3/ Edge Core SWITCH 4598 002 0324x	Figure 1 DACC(2)
	4/ HOS 3T (only 3T cabinets) 4522 132 4457x	Figure 1 DACC(2)
	5/ - Z420 4598 003 2733x / 4598 003 2730x	Figure 1 DACC(2)
	6/ - CDAS ASSY II 4598 000 1447x - CDAS ASSY III 4598 001 8930x	Figure 1 DACC(2)
	7/ - PMU 4535 673 9165x	Figure 1 DACC2
	8 - RF AMPLIFIER 1.5T 4522 132 0324x - RF AMPLIFIER 3.0T 4522 132 0318x - RF AMPLIFIER 3.0T MN 4535 670 7798x	Figure 1 DACC(2)

Batteries 	Items, type	Location
	Not applicable	

Hazardous substances 	Items	Location
	Additional panels added in slots CDAS	Figure 1 DACC(2)
	Panel: RX2 4522 117 9425x	Figure 1 DACC(2)
	Panel: RX3 4522 117 9533x	Figure 1 DACC(2)
	Panel: BULK HUB-2 4522 117 9878x	Figure 1 DACC(2)
	Panel: RX SPLIT. 4535 671 3149x	Figure 1 DACC(2)
	Panel: BULK HUB IF 4535 670 3132	Figure 1 DACC(2)
Panel: BULK HUB 4535 673 2185x	Figure 1 DACC(2)	

Fluids and gases	Items	Location
-------------------------	--------------	-----------------

Name: WEEE RECYCLING PASSPORT <DACC(2)>	12NC: DMR192492	Version: 00 Doc.Date: 2013-12-11
Author:	Page.: 1 of 3	Group No: 5790 Section: - Status: Authorized



Not applicable

Locations as mentioned in the Passport




Figure 1 Frontview
DACC2 cabinet

Identification codes

Name: WEEE RECYCLING PASSPORT <DACC(2)>		12NC: DMR192492	Version: 00 Doc.Date: 2013-12-11
Author:	Page.: 2 of 3	Group No: 5790	Section: -
Status: Authorized			Template version 2.

9896 030 28001	9896 030 27121	9896 030 27071 (LP)	9896 030 26981 (LP)

Name:	WEEE RECYCLING PASSPORT <DACC(2)>	12NC:	DMR192492	Version: 00
Author:		Page.: 3 of 3	Group No: 5790	Doc.Date: 2013-12-11
			Section: -	Status: Authorized