





Product name:	OPERATOR CONSOLE – Z420		
Identification code(s)	9896 030 27131, 9896 030 27091		
Total weight (in Kg)	53 Kg Excluding Computer		
Producer /Manufacturer	Name company:	Philips Medical Systems Nederland B.V.	
	Address:	Veenpluis 4-6	
	Zip code:	5680 DA BEST	
	Country:	The Netherlands	
	Electronic info:		

Recycle info	Substances:	Location
Special attention 	Items:	Location
	Contains PWB materials.	Figure 1 UP TO 4
	Computer Z420 not part of this passport see 4598 004 43331	Figure 1
	Keyboard and mouse's are part of the computer assy.	Figure 2
	AUDIO MODULE 3 (AM3) not part of this passport. See 459800138521	Figure 3
	Display widescreen OCD monitor 24 not part of this passport see 459800069942	Figure 4
Batteries  To be removed	Type:	Location
	NA	
Hazardous  To be removed	Examples	
	NA	
 Fluids / Gases	Items:	Location
	NA	

Name: Recycling Passport Console		12NC: DMR192512		Version: 00 Doc.Date: 2014-01-08	
Author:		Pag.: 1 of 2	Group No: 5790	Section: -	Status: Authorized



Figure 1 Computer Z420



Figure 2 Keyboard and mouse



Figure 3 AUDIO MODULE 3 (AM3)



Figure 4 Display

Name: Recycling Passport Console		12NC: DMR192512	Version: 00 Doc.Date: 2014-01-08	
Author:		Pag.: 2 of 2	Group No: 5790	Section: -
			Status: Authorized	