





Product name	Ingenia 1.5T/3.0T Covers	
Identification code(s)	98960302179.X 98960302180.X 98960302182.X 98960302386.X 98960302619.X	CONES+RINGCOVER REAR WA SIDE+RINGCOVER FRONT WA TOP+TURRET+PLINTH+LOWER COVERS WA MOBILE COVERSET INGENIA PRE-ASSY SYSTEM COVERS WA
Total weight [kg]	168 kg Complete Cover set	
Manufacturer	Philips Medical Systems Nederland B.V.	
Address	Veenpluis 4-6 Best	
Zip code	5680 DA	
Country	The Netherlands	
Electronic info		

Recycling information

Special Attention	Items	Location
	Covers are made of composite materials consisting of glass fibers and resin	Figure 1 and 2
	Rear side: Bore lights (2x) contain both a lamp	Figure 2
	Note: 1.5T and 3.0T have identical cover sets	Figure 1 and 2
	Note: UIMs have their own document	Figure 1

Batteries	Items, type	Location
		

Hazardous substances	Items	Location
		

Fluids and gases	Items	Location
		

Name: Recycling Passport Ingenia 1.5T/3.0T Covers		12NC: DMR176760	Version: 00 Doc.Date: 05-SEP-2013
Author:	Pag.: 1 of 2	Group No: 5790	Section: - Status: Authorized

Locations as mentioned in the Passport (Information in pictures)



Figure 1: Ingenia 1.5T/3.0T front side



Figure 1: Ingenia 1.5T/3.0T rear side

<END OF DOCUMENT>

Name: Recycling Passport Ingenia 1.5T/3.0T Covers		12NC: DMR176760	Version: 00 Doc.Date: 05-SEP-2013
Author:	Pag.: 2 of 2	Group No: 5790	Section: - Status: Authorized