


Product name:	PFEI-WA-iRF 30T		
Identification code(s)	9896-030-1949x		
Total weight (in Kg)	20 Kg		
Producer /Manufacturer	Name company:	Philips Healthcare Best	
	Address:	Veenpluis 4-6	
	Zip code:	5680DA Best	
	Country:	The Netherlands	
	Electronic info:	PFEI-WA-iRF 15T	

Recycle info	Substances:	Location
Special attention 	Items:	
	No special attention needed.	
 Batteries To be removed	Type:	Location
	No batteries present.	
Hazardous  To be removed	Examples	Location
	Depending on the implementation four or six Printed circuit board assemblies	Figure 1a and 1b
	Each Printed circuit board assembly contains 10 times 14mm beryllium copper sticky fingers.	Figure 2
 Fluids / Gases	Items:	Location
	No fluids and/or gases present.	

Name:	WEEE PFEI-WA-iRF 30T, 9896-030-1949	12NC:	DMR115688	Version: 00	Doc.Date: 2010-05-21
Author:		Pag.: 1 of 2	Group No: 5790	Section: -	Status: see Agile

Locations as mentioned in the Passport (Pictures information)




Figure 1a, PFEI-WA-iRF right-side




Figure 1b, PFEI-WA-iRF left-side




Figure 2, Beryllium Sticky-Fingers on Printed circuit board assembly

Note: Figures comprehends the none MN option

Name: WEEE PFEI-WA-iRF 30T, 9896-030-1949		12NC: DMR115688		Version: 00 Doc.Date: 2010-05-21	
Author:		Pag.: 2 of 2	Group No: 5790	Section: -	Status: see Agile
© PHILIPS MEDICAL SYSTEMS NEDERLAND BV, BEST 2010				Template version 2.	