





<b>Product name</b>	<b>DACC-L-D 1.5T</b>
Identification code(s)	See list in below table
Total weight [kg]	250 ± 10 kg
<b>Manufacturer</b>	<b>Philips Medical Systems Nederland B.V.</b>
Address	<i>Veenpluis 4-6 Best</i>
Zip code	<b>5680DA</b>
Country	<b>The Netherlands</b>
Electronic info	WEEE RECYCLING PASSPORT DACC-L-D-1.5T

### Recycling information

<b>Special Attention</b> 	No Special attention required for the following items	<b>Location</b>
	1/ DACC Basic cabinet                      4598 006 5807x	Figure 1 DACC-L-D
	2/MIU-II    4598 006 0652x	Figure 1 DACC-L-D
	6 /TEMP_HUMIDITY_SENSOR    4522 117 9531x	Figure 1 DACC-L-D
	9 /Cooling Interface - Assy            4598 005 9555x	Figure 1 DACC-L-D
	The following items have their own passport and are not part of this passport	<b>Location</b>
	3/ Auxiliary Interface Box (AIBo)    4598 001 3851x	Figure 1 DACC-L-D
	4/ Edge Core SWITCH                      4598 002 0324x	Figure 1 DACC-L-D
	5/ SCC Z230 + DNA-NIC HUB            4598 005 8639x	Figure 1 DACC-L-D
	7/ Transmit control Interface            4598 002 2295x	Figure 1 DACC-L-D
8 /RF AMPLIFIER 1.5T S35            4598 003 6704x	Figure 1 DACC-L-D	

<b>Batteries</b> 	<b>Items, type</b>	<b>Location</b>
	Not applicable	-

<b>Hazardous substances</b> 	<b>Items</b>	<b>Location</b>
	None	-

<b>Fluids and gases</b> 	<b>Items</b>	<b>Location</b>
	Secondary coolant: Distilled or demineralized water with inhibitor and biocide for RF Amplifier circuit	Fluid lines

Name: <b>WEEE RECYCLING PASSPORT &lt;DACC-L-D-1.5T&gt;</b>		12NC: <b>DMR208091</b>		Version: 01 Doc.Date: <b>2014-07-23</b>	
Author:		Page.: <b>1 of 3</b>	Group No: 5790	Section: -	Status: <b>Authorized</b>

<b>Fluids and gases</b>	Properties: Biocide: 0.01 % ( 3.6ml/40l) of NX1164 Inhibitor: 0.02% (8ml/40l) of AZ8104	Fluid lines
-------------------------	---	-------------

**Locations as mentioned in the Passport**




Figure - 1

Name: <b>WEEE RECYCLING PASSPORT &lt;DACC-L-D-1.5T&gt;</b>	12NC: <b>DMR208091</b>	Version: 01 Doc.Date: <b>2014-07-23</b>
Author:	Page.: <b>2 of 3</b>	Group No: 5790
Section: -	Status: <b>Authorized</b>	

## Identification codes

98960302708x	DACC-L-D-1.5T
--------------	---------------

Name: <b>WEEE RECYCLING PASSPORT &lt;DACC-L-D-1.5T&gt;</b>			12NC: <b>DMR208091</b>			Version: 01 Doc.Date: <b>2014-07-23</b>		
Author:			Page.: <b>3 of 3</b>	Group No: <b>579O</b>	Section: <b>-</b>		Status: <b>Authorized</b>	