



Alice NightOne home sleep testing system

Alice NightOne advantages

- “Auto-on” feature
- RIP with integrated buckle design
- Unique “smart guide” patient setup sequence
- Enhanced Good Study Indicator
- Body Position sensor built-in
- Powered by Sleepware G3 software
- Optional Cloud services
- Signals your techs are familiar with
- Wireless connection to System One PAP device

Designed for ease, confidence and reliability

The Alice NightOne home sleep testing (HST) solution combines ease, comfort and reliability to support your business, clinical and patient needs when it comes to HST. With features and functions designed with patients in mind, Alice NightOne is up to two times more likely to result in an acceptable study than competitive home sleep testing devices.*

Alice NightOne includes:

- Alice NightOne device
- Hard shell carrying case
- Nasal cannula sample pack
- zRIP effort belt
- Oximetry module
- SpO₂ sensor (style may vary)
- USB cable
- Quick Start guide, user and provider manuals
- Sleepware G3 software
- Unlimited software updates and site licensing
- Two-year device warranty



*Results of customer preference trial; data on file.

Alice NightOne information

Specifications

Total channels	Seven, plus optional therapy channels
Recorded channels	Flow, snore, thoracic effort, SpO ₂ , pulse rate, body position, pleth; plus optional therapy channels depending on therapy device and mode
Body position reporting	Supine, prone, right side, left side, upright
Good Study Indicator reporting	Displays quality and quantity of signals recorded, allowing clinicians to identify/correct poor signals
Data storage	4 gigabyte internal memory
Communication	USB to PC cable, wireless communication to Bluetooth-enabled Philips Respironics devices
Power requirements	Two AA (1.5V) alkaline batteries or two AA (1.2V) rechargeable batteries with at least 2400mAh capacity
Oximetry sampling rate	62.5 Hz
Dimensions	4.07" L X 0.99" W X 2.67" H (10.34 cm X 2.51 cm X 6.78 cm)
Weight	Approximately 3 oz. (84 grams); not including sensors or batteries
Easy to clean	Choice of mild soapy water, Isopropyl alcohol (70%), DisCide®, CaviCide® or EnviroCide®
Device warranty	Two years

Part numbers

Device	Part number	Replacement accessories	Part number
Alice NightOne device kit (US/Canada)	1113277	Thoracic effort belt, adult	1113310
		Thoracic effort belt, large	1113311
		Pro-Flow nasal cannula, adult, qty. 60	P1328-60
		Carrying case	1117455
		Case insert replacement kit	1114520
		Quick Start guide, 10 pack	1113329

