



**PHILIPS**

*DreamStation Go*

Portable PAP therapy system

# Innovation

for people on the Go



© 2018 Koninklijke Philips N.V. All rights reserved.  
Specifications are subject to change without notice.

[www.philips.com.sg/sleepapnea](http://www.philips.com.sg/sleepapnea)

DC 10/9/2018 PN1867143

# Convenience without compromise

The **DreamStation Go** is a compact and highly portable PAP system designed for patients who want to maintain their therapy at home and on the road.



#### Powerful, portable performance

At half the size of our System One and DreamStation platforms, DreamStation Go offers the same therapy and comfort features used by more than 5 million CPAP users.



#### Mask compatibility

DreamStation Go gives your patients the freedom to choose the mask they want.<sup>1</sup>



#### Optional battery

DreamStation Go's overnight battery can help your patients stay committed to their therapy when they are off the grid. Set at 10cm of pressure, DreamStation Go displayed an average run time of 13 hours<sup>2</sup> when operating on its overnight battery.



#### 12mm micro-flexible tubing

DreamStation Go features the smallest PAP tubing we've ever offered, providing enhanced flexibility, portability, and connectivity. In a recent head-to-head trial, patients found the DreamStation 12mm micro-flexible tubing to offer better flexibility, more comfort, and greater ease of movement in comparison to standard tubing.<sup>3</sup>



#### Simple setup

With our intuitive color touch screen navigation, a mobile device app isn't required to set up or adjust your patients' therapy.



#### Compliance reporting made easy

Like our other PAP devices, DreamStation Go is compatible with our patient data management system, which seamlessly combines usage and compliance data from multiple Philips Respironics CPAP devices into one simple report.



#### DreamMapper cross-device connectivity

Users of DreamMapper, our patient self-management tool, can sync data from multiple DreamMapper-enabled devices<sup>4</sup> to one account. They also have the ability to track and monitor their combined progress across both home and portable PAP devices via one convenient dashboard on their smartphone, tablet, or computer.

<sup>1</sup>Masks with standard 22mm outer diameter connection.

<sup>2</sup>Battery run times may vary due to higher CPAP pressure settings, excessive mask leakage, altitude, temperature, charging devices with the USB port, and other variables.

<sup>3</sup>Head-to-head trial of DreamStation 15mm tubing vs. DreamStation 12mm micro-flexible tubing, 2016.

<sup>4</sup>To see which therapy devices are compatible with DreamMapper, visit [www.dreammapper.com/compatible](http://www.dreammapper.com/compatible).





# The confident, convenient choice

DreamStation Go provides the portability and packability travelers want and the effective therapy they need.

**65%  
smaller\***

The DreamStation Go is 65% smaller than ResMed's AirSense 10 and AirStart 10.

**58%  
smaller**

The DreamStation Go 12mm tubing is 58% smaller than our 22mm performance tubing.

**48%  
lighter\***

The DreamStation Go is 48% lighter than ResMed's AirSense 10 and AirStart 10.

**33%  
fewer components\*\***

The DreamStation Go requires 33% fewer components to pack than its closest competitors.

\*Internal measurement of ResMed's AirSense 10 and AirStart 10 with side cover installed. Includes power supply. Data on file, 2016.

\*\*When compared to HDM Z1 and Sommetics Transcend devices, 2016.

# Powerfully convenient features



Optional battery provides approximately 13 hours of therapy!<sup>1</sup>



Intuitive color touchscreen offers quick, simple setup and easy navigation.



Built-in power supply reduces the number of parts managed by the patient.



Optional microSD card stores compliance and waveform data.<sup>2</sup>



DreamStation Go is the only travel PAP with a USB charging port,<sup>3</sup> allowing for easy charging of devices.

<sup>1</sup>PAP run-time test conditions: CPAP mode, pressure 10cm H<sub>2</sub>O, 12mm tubing, 37LPM of leak, 73.4 F room at an elevation of 1650 ft.

<sup>2</sup>Sold separately.

<sup>3</sup>When compared against HDM Z1, Somnatics Transcend, and ResMed's AirMini devices, 2017.

<sup>4</sup>RTCA/DO-160G section 21, category M.



Bluetooth is standard with every device.



Unit complies with FAA requirements for in-flight use.<sup>4</sup>



Device is compatible with DreamMapper, our patient self-management tool.

## 12mm tubing

Micro-flexible tubing is our smallest and lightest ever.

## Click-connections

Click-connections on each end of the micro-flexible tubing make connecting and disconnecting the device and mask simple.





# Create a world of opportunities



**Overnight battery**  
Compact and lightweight, the battery features a recharge time of just 4-5 hours.



**Small travel kit**  
Kit includes an easy-to-pack protective case for the DreamStation Go and a washable bag for the mask and tube.



**Medium travel kit**  
Kit includes a protective case that accommodates the DreamStation Go and overnight battery, as well as a washable bag for the mask and tube.



**Air filters**  
Choice of a reusable pollen or disposable 30-day fine filter.



**12mm micro-flexible tubing**  
This tubing offers more comfort and greater ease of movement than that of standard tubing.\*

\*Head-to-head trial of DreamStation 15mm tubing vs. DreamStation 12mm micro-flexible tubing, 2016.

## DreamStation Go Auto specifications

Pressure range	4 to 20 cm H <sub>2</sub> O
Pressure mode(s)	CPAP, APAP
Flex pressure relief	1 to 3
Flex modes	C-Flex, C-Flex+, A-Flex
Ramp time	0 to 45 min (5-min Increments)
Starting ramp pressure	Patient adjustable, Smart Ramp
Dimensions	15 cm L x 15 cm W x 5.8 cm H
Weight	844 g (1.86 lb)
Filters	Reusable pollen, combination fine
Device setup	Touchscreen color display
Data storage capacity	On-board: > 1 year MicroSD card (optional): lifetime
Compliance meter	Breathing detection
Altitude compensation	Automatic
Electrical requirements	100-240 VAC, 50/60 Hz
Warranty	Two years

## Ordering information

Flow generators	Part number
DreamStation Go Auto CPAP	Contact your local sales office
Power	Part number
DreamStation Go overnight battery	1133281
Travel kits	Part number
DreamStation Go small travel kit*	1133275
DreamStation Go medium travel kit**	1133276
Consumables/Replacement	Part number
DreamStation micro-flexible 12mm tubing	PR12
Power cord 6ft white	1132551
Power cord 10ft white	1132552
Beauty panel - front	1133741
Beauty panel - rear	1133742
Reusable pollen filter 2-pk	1133743
Disposable 30-day fine filter 6-pk	1133744



## Learn more

For more information about DreamStation Go, visit [www.philips.com/dreamstationgo](http://www.philips.com/dreamstationgo).

Toll free number

Philippines 1800-1888-6182 | Singapore 1800-28-63-020 | Thailand 1800-999-119

\*Small travel kit includes an easy-to-pack protective case for the DreamStation Go and a separate washable mask and tube bag.

\*\*Medium travel kit includes a protective case that accommodates the DreamStation Go and overnight battery. Also includes a separate washable mask and tube bag.