





Product name:	PFEI-C-10T		
Identification code(s)	4522-132-2023		
Total weight (in Kg)	20 Kg		
Producer /Manufacturer	Name company:	Philips Medical Systems Best	
	Address:	Veenpluis 4-6	
	Zip code:	5680DA Best	
	Country:	The Netherlands	
	Electronic info:	PFEI-C10T_452213220230_579_v3.doc	

Recycle info	Substances:	Location
Special attention 	Items:	
	No special attention needed.	
Batteries  To be removed	Type:	Location
	No batteries present.	
Hazardous  To be removed	Examples	Location
	Depending on the implementation four or six Printed circuit board assemblies	Figure 1a and 1b
	Each Printed circuit board assembly contains 10 times 14mm beryllium copper sticky fingers.	Figure 2
 Fluids / Gases	Items:	Location
	No fluids and/or gases present.	

Name: WEEE passport, PFEI-C-10T		12NC: 4522-132-2023		Version: 0 Doc.Date: 27-07-2006	
Author: R.Kunnen		Pag.: 1 of 3	Group No: 579	Section: 0	Status: CONCEPT

Locations as mentioned in the Passport (Pictures information)




Figure 1a, PFEI-C-10T right-side




Figure 1b, PFEI-C-30T left-side




Figure 2, Beryllium Sticky-Fingers on Printed circuit board assembly Recycling Passport Revision History

Name: WEEE passport, PFEI-C-10T		12NC: 4522-132-2023	Version: 0 Doc.Date: 27-07-2006
Author: R.Kunnen	Pag.: 2 of 3	Group No: 579	Section: 0
			Status: CONCEPT

Recycling Passport Revision History

Rev	Description	Date
1	First Passport (Draft)	March 31, 2006
2	Second Version (Concept)	27 July, 2006
3	Third Version (Authorised)	28 July, 2006

Recycling Passport Approval Signatures

Written By R.Kunnen

R.Kunnen

Reviewed By (optional)

M.Boutelje

Reviewed By (optional)

D.Curiel

Approved By

GJ Plattel

Name: WEEE passport, PFEI-C-10T		12NC: 4522-132-2023	Version: 0 Doc.Date: 27-07-2006
Author: R.Kunnen	Pag.: 3 of 3	Group No: 579	Section: 0 Status: CONCEPT