





Product name:	Allura XPER FD20		
Identification code(s)	4522 901 04440		
Total weight (in Kg)	2000 Kg (approx)		
Producer /Manufacturer	Name company:	Philips Medical Systems	
	Address:	Veenpluis 6	
	Zip code:	5684 PC Best	
	Country:	The Netherlands	
	Electronic info:	Allura_FD20_579.doc	

Recycle info	Substances:	Location
Special attention 	Items:	Location
	Be aware of possible contaminated system parts and materials.	See Figure 1
	For dismantling activities Treatment Facilities must consider the national requirements.	
	For personnel that can come in contact with contaminated material, preventive measurements pursuant to national requirements must taken into account.	
	Removal of units/weights can cause the system to tilt!	
	Removal of units/weights can cause unexpected movements.	
	High voltage parts are marked with	
Hazardous 	Lead (Pb) for X-Ray Shielding	
	Lead (Pb) for soldering	All electronics
	Beryllium copper (BeCu)	Some electronics
	Mercury (Hg) in specific LCD Screens	LCD Screens
To be removed		
 Fluids / Gases	Items:	Location
	Transformer Oil, type	

Note: To facilitate recycling, all plastic parts > 50 grams are marked according to ISO 11469 & ISO 1043

Name:	WEEE Passport Allura XPER FD20	12NC:	4522 901 04440	Version: 00
Author: W.Bleyen		Pag.: 1 of 3	Group No: 579	Doc.Date: 2006-09-01
			Section: -	Status: Approved

Figure 1



See Frontal Stand
R-Passport

See MCS
R-Passport

See Detector
R-Passport

See User Interface
R-Passport

See Patient
Table R-Passport

See X-Ray Tubes R-Passport
See Collimator R-Passport
See Grid R-Passport

For Stand:
See Frontal Stand Allura XPER FD20

For Cabinets :
See R-Cabinet R-Passport
See M-Cabinet R-Passport
See X-Ray Generator Passport for Velara CV Cabinet

Name: WEEE Passport Allura XPER FD20		12NC: 4522 901 04440		Version: 00 Doc.Date: 2006-09-01	
Author: W.Bleyen		Pag.: 2 of 3	Group No: 579	Section: -	Status: Approved

Recycling Passport Revision History

Rev	Description	Date
1	First Passport	Sept. 1, 2006

Recycling Passport Approval Signatures

Written By <b style="color: blue;">PMS W.Bleyen
--

Reviewed By (optional)

Approved By <b style="color: blue;">PMS
--

Approved By

Name: WEEE Passport Allura XPER FD20	12NC: 4522 901 04440	Version: 00 Doc.Date: 2006-09-01
Author: W.Bleyen	Pag.: 3 of 3	Group No: 579
	Section: -	Status: Approved