





Product name:	PX4800	
Identification code(s)	9896 010 23491, 9896 001 85001.2.3.4, 9896001, 9896 001 86331	
Total weight (in Kg)	10.3 +/- 0.2 kg	
Producer/ Manufacturer	Name company:	Trixell S.A.
	Address:	460, Rue du Pommarin
	Zip code:	38 430, Moirans
	Country:	France

Recycle Info	Items:	Location
Special attention 		
Fluids / Gases 	Items:	Location
	Glycoshell SF, 58 ml (pink fluid)	In plumbing (cooling plate and tubes)
Batteries 	Type:	Location
Hazardous  To be Removed	Substances:	Location
	Lead (Pb), 2995 grams	See figure
	Thallium, max 1,9 grams, used in the Scintillator that converts Xrays into light.	See figure

Printed copies are uncontrolled

Title: _____
 Recycling Passport Number: _____ Rev: _____

Show locations of materials mentioned on the previous sheet. (Picture information)

Printed copies are uncontrolled




From : I.Quesse

Date : 06/06/06

Ref. : Tri 06/06/0705

COMMUNICATION

Lead and Thallium parts in Pixium 4800



Copyright © 2006 TRIXELL. Tous droits réservés. Ce document et l'information qu'il contient sont la propriété exclusive de TRIXELL. Il ne peut être reproduit ou divulgué sous quelque forme que ce soit sans l'accord préalable écrit de TRIXELL. TRIXELL PROPRIETARY INFORMATION. Copyright © 2006 by TRIXELL. All rights reserved. This document and the information contained therein are TRIXELL's exclusive property. Any copy and/or disclosure thereof in any form whatsoever are subject to TRIXELL's prior written consent.

The plate just right to the utmost right lead plate is the cooling plate. This contains 58 ml Glycoshell SF.

Title: _____

Recycling Passport Number: _____ Rev: _____

Printed copies are uncontrolled

Recycling Passport Revision History

Rev	Description	Date

Recycling Passport Approval Signatures

<p>Written By</p> <p>Marcel van der Steen</p>
--

<p>Reviewed By (optional)</p>

<p>Approved By</p> <p>Joost Jacobs</p>

<p>Approved By</p>

Title: _____

Recycling Passport Number: _____ Rev: _____