

Product name:	Patient Table FD10_ FD20		
Identification code(s)	4522 901 00510		
Total weight (in Kg)	400 Kg (approx)		
Producer /Manufacturer	Name company:	Philips Medical Systems	
	Address:	Veenpluis 6	
	Zip code:	5684 PC Best	
	Country:	The Netherlands	
	Electronic info:	Patient Table FD10_FD20_579.doc	

Recycle info	Substances:	
Special attention 	Items:	Location
	Be aware of possible contaminated system parts and materials.	See Figure 1
	For dismantling activities Treatment Facilities must consider the national requirements.	
	For personnel that can come in contact with contaminated material, preventive measurements pursuant to national requirements must taken into account.	
	Removal of units/weights can cause unexpected non motorized movements!	
	Removal of units/weights can cause the system to tilt!	
	High voltage parts are marked with	
Hazardous 		
	Lead (Pb) for soldering	All electronics
	Beryllium copper (BeCu)	Some electronics
To be removed		

Note: To facilitate recycling, all plastic parts > 50 grams are marked according to ISO 11469 & ISO 1043


Name:	WEEE Patient Table FD10_ FD20	12NC:	4522 901 00510	Version: 00
Author: W.Bleyen		Pag.: 1 of 3	Group No: 579	Doc.Date: 2006-09-01
			Section: -	Status: Approved

Patient Table

Figure 1



Figure 2



Beryllium copper (BeCu)

Recycling Passport Revision History

Rev	Description	Date
1	First Passport	Sept. 1, 2006

Recycling Passport Approval Signatures

Name: WEEE Patient Table FD10_FD20		12NC: 4522 901 00510		Version: 00 Doc.Date: 2006-09-01	
Author: W.Bleyen	Pag.: 2 of 3	Group No: 579	Section: -	Status: Approved	

Written By
PMS K.A.J. de Graaff

Reviewed By (optional)

Approved By
PMS

Approved By

Name:	WEEE Patient Table FD10_FD20	12NC:	4522 901 00510	Version:	00
Author:	W.Bleyen	Pag.:	3 of 3	Doc.Date:	2006-09-01
		Group No:	579	Section:	-
		Status:	Approved		