


Printed copies are uncontrolled


Product name:	M-Cabinet FD	
Identification code	9896 010 4031* : M-CABINET FDXD 9896 010 4044* : M-CABINET FDXD VASCULAR 9896 010 4048* : M-CABINET FD MONOPL 9896 010 4049* : M-CABINET FD EXT. TO BIPLANE 9896 010 4056* : UPGRADE PBL10 TO PBL30	
Total weight (in kg)	300-340 kg	
Producer/ Manufacturer	Name company:	Philips Medical Systems B.V.
	Address:	Veenpluis 4-6 Best
	Zip code:	5508AD
	Country:	The Netherlands
	Electronic info:	

Recycle Info	Items:	Location
Special attention 		
Fluids / Gases 		
Batteries  To be Removed	Type:	Location
	Alkaline battery	1
Hazardous  To be Removed	Substances:	Location
	BeCu Contact fingers	2 Dotted area (at FRONT- and BACKSIDE cabinet)

Title: Eco_passp_M-cabinet
 Recycling Passport Number: 8222 168 61590 Rev00

Show locations of materials mentioned on the previous sheet.

Printed copies are uncontrolled



Recycling Passport Revision History

Rev	Description	Date
00	New	June 23 2006

Recycling Passport Approval Signatures

Written By
W van den Berg

Reviewed By

Approved By
B. Gillis

Approved By

Printed copies are uncontrolled